

equipment &
consumables

TASK
MEDICAL



taskmedical.com.au



X ray gown

- Belt loop
- Disposable
- Shoulder sleeve
- Generous size
- Box of 100

Standard X large
GOWNT GOWNXL

Isolation gown

- Belt loop
- Long sleeve
- Disposable
- Box of 100

Standard
GOWNY



Face mask

- High filtration
- 3 ply
- Ear loops
- Pack of 50

TK0301



Task Medical Tourniquet Quality product made in Germany

Proven for daily use in hospitals, medical practices, and for paramedics.

- Skin-friendly, flexible and soft ribbon 25 mm wide
- Sturdy safety snap closure
- One-hand operation
- Metalfree
- No squeezing or pinching of the skin
- Latex free

TKTOUR



Head set covers

- Elastic rim
- Disposable
- Pack of 250
- Hygienic

HEA



Fitted couch cover



- Non woven
- Disposable
- 200cm x 75cm
- Carton of 100

White
TKFITWHT
 Navy Blue
TKFITSNB


Stretcher sheet



- Non woven
- 240cm x 70cm
- Disposable
- Carton of 100

White
SHEETWHT
 Navy Blue
SHEET

Pillow sleeves



- Disposable
- 70cm x 45cm
- Disposable
- Carton of 200

White
PILWHT
 Navy Blue
PILBLU

TASK MEDICAL PREMIUM

sms 40gsm
 non woven fabric

Task Medical premium range consumables are made from a superior 40gsm SMS non woven fabric, comprising of two spunbond layers interlaid with one meltblown layer. They are an effective hygienic barrier against contamination on examination tables protecting your patients and furniture surfaces. Tear resistant and breathable.

Fitted bed sheets



- 200cm x 70cm x 10 cm
- Disposable
- Fitted elastic ends
- Fitted corners
- Fluid resistant

Light blue
TKPRFITS-LB

Modesty sheets



- 100cm x 70cm
- Disposable
- Patient cover
- General purpose
- Fluid resistant

Light blue
TKPRMODS-LB

Pillow cases



- 70cm x 45 cm
- Disposable
- Sealed on 3 sides
- Open end flap
- Fluid resistant

Light blue
TKPRPLCS-LB

Versatile, popular product



TASK Multi Roll Available in 3 sizes

- Highly absorbent
- 2ply Perforated
- 25cm x 50m
- 12 rolls per carton

- Highly absorbent
- 2ply Perforated
- 50cm x 50m
- 6 rolls per carton

- Highly absorbent
- 2ply Perforated
- 58cm x 50m
- 6 rolls per carton




TKMR25

TKMR50

TKMR58

TASK Bed Roll


- Keep furniture & surfaces clean between patients
- Added protection
- 1 ply
- 55cm x 80m
- 6 rolls per carton



TKBR55

TASK Utility Roll - Laminated

- Fluid proof & hygienic barrier
- Perforated Polyethylene Backed
- 68cm x 50m
- 12 rolls per carton



TKUR6850

Clarity Ultrasound Gel

- Not affected by the salt on the skin during the usage
- Ensures a comfortable ultrasonography
- Hypoallergenic, does not irritate the skin
- Doesn't contain oil or oily agents, formaldehyde or salt
- No toxic effect or odour

ALSO AVAILABLE

Clarity Sterile Ultrasound Gel 20ml

Sterilised with gamma rays, this gel is used for ultrasonography and biopsy sampling processes in vitro fertilization, surgeries and vaginal and intra uterine examinations.



TASK Medical Probe Covers For Ultrasound

Specially designed for Gynaecological Ultrasonography Each probe cover is rolled for easy application to cover a variety of trans-vaginal transducers. 100 individually foil sealed probe covers in a dispenser box.

- Transducer probe cover specifications
- High quality homogenised low allergy rubber latex
- Meets Quality Standard ISO 4070
- Single use / Non Sterile
- L: 200mm
- W: 53mm
- Dia: 34mm
- Available in Lubricated & non-Lubricated

Clarity Probe Cleaning Solution 250ml

Ultrasceptic ultrasound probe cleaning spray solution removes the gel left over the probe. Helps prevents contamination from patient to patient thanks to its active agents and ensures hygiene.



TASK Modesty sheet

- 120cm x 80cm
- 3ply Disposable
- Case of 100
- Protect & Respect



XL Size
TKMS80120

Basic dressing pack

- 1 x Sterile field
- 1 x Dressing tray
- 1 x Dressing towel
- 5 x Cotton balls
- 2 x Forceps blue
- 1 x Forceps yellow
- 2 x Gauze swabs 7.5cm x 7.5cm



TKBDP **QUALITY** Tested *****
STERILE EO 0197

Underpads


TASK underpads are ideal for fluid management, providing extra protection for beds, trolleys, trays etc.

- High wicking
- Fibrous pulp
- Leak prevention
- Highly absorbent
- Sealed 4 sides

Premium
60cm x 40cm
Ctn 300
TKUP30

Economy
60cm x 40cm
Ctn 300
TKUP22

Half Size
Economy
30cm x 40cm
Ctn 400
5 ply
TKUP25



fibrous pulp is
20% more
absorbent
equivalent to 8 ply

Everyday wipes

- Highly absorbent
- Dry wipes
- Latex free
- 30cm x 33cm
- 50 wipes / pack
- Carton of 12 packs



TKEW50

Head pads

- Hygienic barrier
- Pillow protection
- 50cm x 30cm
- 3 ply
- Carton of 1000



TKHP50

Relieve back stress



- Improved posture
- Improved comfort
- Better ergonomics
- Premium Brand
- Lower Prices



Mica Saddle Seat

Tilting seat, easy to change height position by using handle

WINBEX
World Class Quality

Authorised Australian Distributor

Navy Blue
WX5113201NB

Black
WX113201BL

Premium quality stools at affordable prices



Premium Round Stool

Single lever for height adjustment

Navy Blue
WX002327NB
Black
WX002327BL
Grey Teal
WX002327GT



Comfort Line Stool

140mm / 200mm lift adjustment

Navy Blue
WX002908NB
Black
WX002908BL
Grey Teal
WX002908GT



Ascot Standard Saddle Seat

Design and built for ergonomics and comfort


Navy Blue
WX5105S02NB
Black
WX5105S02BL



Kimi Saddle Seat

Helps to maintain correct posture

Navy Blue
WX511201NB
Black
WX511201BL



Doctors Stool T Seat

Comfort seat and backrest supports body and spine

Black
WX510710BL



Duma Doctor's Seat

Fully supported body and spine in a sedentary position


Navy Blue
WXDUMANB



Ascot Standard Saddle Seat Back

Backrest provides more support for longer working times


Navy Blue
WX511114NB
Black
WX511114BL



Kimi Saddle Seat with Back

Relieve stress from knees and hips

Navy Blue
WX511202NB
Black
WX511202BL



Microscope Chair

Separate adjustment for seat height, seat tilt and back angle

Black
WX5113301BL



e200 transport stretcher

Equipped with a four-wheel brake system and a central fifth wheel which offers excellent mobility to provide the best safety for patients, the e200 provides durable quality and a 300kg weight capacity.

TKST200



- Foot pedal for hydraulic system control on both sides
- Integrated storage tray
- Two-section bed board / Collapsible side rails
- Drainage Bag clasp x 2 / Urinary Bag clasp x 2
- Foot controlled brake / Retractable steering fifth wheel
- Steer pedal for four wheels at both ends
- Head-end delivery push bar / Foot-end assist push bar
- Dual pneumatic assisted backrest
- Two stage movable IV pole x 1
- Two stage permanent folding IV pole x 1
- Velcro tape underneath the mattress
- Protective corner bumpers x 4
- Hook / jack for upright oxygen tank holder
- 300 kg SWL

Multi Purpose Mobility Bed

EZYM HI LO
by TASK MEDICAL

- Recovery room
- Procedure room
- Patient transporter
- Infusion therapy bed
- Emergency stretcher
- 180kg SWL
- Manual wind
- 56-86cm

Ideal for any medical practice!

TKEZYM



Easy roller mri safe wheel chair

MR SAFE Certified MR safe (7 Tesla) by Siemens and Phillips

Made from recyclable PE and PA thermoplastics

100% non metallic non magnetic

EasyRoller has been especially designed as a MR SAFE patient transport device for use in hospitals and MR laboratories. Patients can easily be moved through the entire hospital and all the way through the MRI room without interfering with the MRI equipment. Norwegian design, German made.



- No metal parts - Non metallic
- No maintenance recyclable plastic
- Easy maneuvering swivel wheel axle
- Tubeless soft ride tyres & park brake
- Double POM glass ball bearing wheel axles
- Waterproof padded seat, back & armrests
- Approved Standard EN12183 and ISO 7176



Transport patients directly to the mri machine

MRI Safe Wheelchair

WHEEASYRL



Single no drawer



TK0D50N

Single one drawer



TK1D50N

Single two drawer



TK2D50N

Single three drawer



TK3D50N

Single four drawer



TK4D50N

Mayo instrument



TKMT3L

Stainless steel trolleys

Task Medical trolleys are made from high grade 304 stainless steel throughout, with double safety folds and smooth welding.

The top and bottom shelves are surrounded by safety rails on all three sides.

Movability is made easy by 95mm smooth rolling castors with lock brakes on the front.

The drawers are soft closing magnetic with German made ball bearing runners.



TK0D80N



TK2D80N



TK4D80N

Emergency Trolley

Have your CPR & Resus' gear ready to go in the one unit

TKEM65032



Clinical Trolley

TKCT1D

Extra large castors for effortless maneuverability



Camino Medical Hi-Lo **TASK** MEDICAL


Electric Exam Couch With/without 360° walk-free foot bar control 3 section - Single Motor

This is a Medical Exam Couch with exceptional strength and stability which outperforms many in its class

- Foot pedal control
- Face hole
- Thicker 64mm foam padding
- Gas strut tilt control
- Steel frame
- Front bar key lock
- Lower & raise retractable castors
- Adjustable feet
- Fireproof PU/leather
- 250kg SWL

- Black
102531
- Charcoal
102531
- Agate blue
102360

- Walk free control bar
Agate blue
102361



Superior German made motor

250kg SWL

WARRANTY
5 YEAR FRAME
3 YEAR ELECTRICAL



Camino Physio Hi-Lo Electric Treatment Table

Designed to support multiple posture adjustments for physio treatments

- Electric height adjust.
- 46-95cm. SWL 250kg
- A double-sided synchronous pneumatic control of the tail-adjustment mechanism makes it more sturdy, and achieves quiet and fast adjustment

Dual struts



Lower & raise retractable castors

Longer section

Agate blue
102535




Sonora All Electric Medical Aesthetic Chair

Beauty & Cosmetic 3 section - 3 Motor Available in Black or White

The Sonora contributes to serve beauty customers with more comfort and indulgence and to improve the treatment capabilities and convenience of beauticians and medical aesthetic clinicians.

- Electric height adjust.
- 56-90cm. SWL 200kg
- Electric backrest
0°-75° angle adjustment
- Electric leg rest
0°-90° angle adjustment
- Body Oil resistant PU leather
- Luxurious 82mm thick high density
- CA-117 fireproof foam padding
- Removable leveled side armrests
- Removable /adjust' neck/face pillow



3 Motors
Height motor
Backrest motor
Footrest motor

Height & multi-position hand controller

Stylish mono frame lift mechanism

Sturdy base

- Black
102363
- White
102362

Face hole

180° swing back arms

Superior German made Motors

More positions

3 Motor control

Adjustable foot rest extension



Dental Bib 2+1 & 1+1

All colours and patterns are available in 2+1 and 1+1



Saliva Ejectors



Head Rest Cover



Dental Roll 30cm x 50m



Dental Cups 180ml

Now available in a range of colours



Dental Cotton Rolls 300gr



Medicine measures

30ml

Clear

MM30

Quantity
5000 pieces/carton;
(50 sleeves x 100 pieces)

Size
30ml - 45mm x 34mm (approx.)



60ml



Clear
MM60

Quantity
1000 pieces/carton;
(40 sleeves x 25 pieces)

Size
60ml - 45mm x 55mm (approx.)

Medicine measure lid
Clear
MML30

Gallipots

60ml

Quantity
1000 pieces/carton

Size
60ml - 70mm x 25mm (approx.)

Clear

BG60C



150ml

Quantity
1000 pieces/carton

Size
150ml - 77mm x 55mm (approx.)



Clear
BG150C

Denture Bowls & Lids

250ml Denture bowl & lid

Quantity
50 pieces/carton

Size
250ml - 95mm x 53mm (approx.)



Clear
B250CL50
Blue
B250BL50
Green
B250GL50

Kidney dishes

700ml single use

Quantity
250 pieces/carton;
(10 bags x 25 pieces)

Size
700ml - 230mm x 120mm x 45mm (approx.)



Clear

DK700CS

Blue

DK700BS

Green

DK700GS

Yellow

DK700YS

700ml autoclavable

Quantity
100 pieces/carton;
(10 bags x 10 pieces)

Size
700ml - 230mm x 120mm x 45mm (approx.)



Clear

DK700CU

Blue

DK700BU

Green

DK700GU

Yellow

DK700YU

1000ml autoclavable

Quantity
100 pieces/carton;
(10 bags x 10 pieces)

Size
1000ml - 320mm x 150mm x 50mm (approx.)



Clear

DK1000CU

Blue

DK1000BU

Green

DK1000GU

Yellow

DK1000YU

Disposable procedure trays

200ml

Quantity
350 pieces/carton;
(7 bags x 50 pieces)

Size
200ml - 200mm x 70mm x 30mm (approx.)



Clear

TP200C

Blue

TP200B

Green

TP200G

Yellow

TP200Y

500ml

Quantity
175 pieces/carton;
(7 bags x 25 pieces)

Size
500ml - 260mm x 90mm x 40mm (approx.)



Clear

TP500C

Blue

TP500B

Green

TP500G

Yellow

TP500Y

1500ml

Quantity
80 pieces/carton;
(4 bags x 20 pieces)

Size
1500ml - 245mm x 205mm x 40mm (approx.)



Clear

MT01C

Blue

MT01B

Green

MT01G

Yellow

MT01Y

Disposable bowls

250ml

Quantity
500 pieces/carton;
(20 bags x 25 pieces)

Size
250ml - 95mm x 50mm (approx.)

Clear Blue Green
B250C B250B B250G



500ml

Quantity
300 pieces/carton;
(12 bags x 25 pieces)

Size
500ml - 115mm x 60mm (approx.)

Clear Blue Green
B500C B500B B500G



1000ml

Quantity
270 pieces/carton;
(9 bags x 30 pieces)

Size
1000ml - 145mm x 80mm (approx.)

Clear Blue Green
B1000C B1000B B1000G



2000ml

Quantity
100 pieces/carton;
(4 bags x 25 pieces)

Size
2000ml - 195mm x 90mm (approx.)

Clear Blue Green
B2000C B2000B B2000G



5000ml

Quantity
50 pieces/carton;
(5 bags x 10 pieces)

Size
5000ml - 260mm x 110mm (approx.)

Clear Blue Green
B5000C B5000B B5000G



6000ml

Quantity
40 pieces/carton;
(4 bags x 10 pieces)

Size
6000ml - 295mm x 110mm (approx.)

Clear Blue Green
B6000C B6000B B6000G



Autoclavable bowls

150ml

Quantity
10 pieces/carton

Size
150ml - 100mm x 50mm (approx.)

Yellow Green
B100YU B100GU



500ml

Quantity
10 pieces/carton

Size
500ml - 140mm x 60mm (approx.)

Green Yellow
B140GU B140YU



1000ml

Quantity
10 pieces/carton

Size
1000ml - 185mm x 90mm (approx.)

Green Yellow
B185GU B185YU



1500ml

Quantity
10 pieces/carton

Size
1500ml - 205mm x 80mm (approx.)

Yellow Green
B205YU B205GU



2500ml

Quantity
10 pieces/carton

Size
2500ml - 240mm x 110mm (approx.)

Green Yellow
B240GU B240YU



3500ml

Quantity
10 pieces/carton

Size
3500ml - 305mm x 90mm (approx.)

Yellow Green
B305YU B305GU



6000ml

Quantity
10 pieces/carton

Size
6000ml - 345mm x 120mm (approx.)

Yellow Green
B345YU B345GU

