



Clinical Range

TASK
MEDICAL

***All your basics from
a single source***

Non-Woven Bouffant Cap

White
TKBFCPW

- White/Blue
- PPSB 18g/m²
- 21" (A:30cm+/-2cm)
- Round Elastic, Latex Free

Blue
TKBFCPB



Non-Woven Mob Cap

White
TKMBCPW

- White/Blue
- PP 12gsm
- 21"
- Double Elastic, Latex Free

Blue
TKMBCPB



Non-Woven Beard Cover

- Blue
- PPSB12-14gsm
- 18" Latex Free
- Nylon, Polyester

Blue
TKBDCV



Non-Woven Shoe Cover

- Blue
- Anti-Skid
- 17 X 40cm
- Textured for better grip

Blue
TKSHCV



Isolation Gown

Blue
TKGOWNISOB

- Light Blue
- PP+PE 30g+15
- Disposable
- Splash Resistant
- Micro Breathable M



Knitted Cuffs

Lab Coat

- White
- PPSB 40gsm
- 120cm x 145cm
- Knitted cuffs
- 4 Velcro fasteners

White
TKLBCT

Single Collar



W/O Pockets



Examination Gown Sleeveless

- Patient comfort
- Patient modesty
- Dark Blue - Opaque
- 2 sizes available
- Ties

Order Code
GOWNT

Order Code
GOWNXL



SMS Isolation Gown

Order Code
TKISGNSMSB

- Light blue
- 120 X 140cm
- Knitted Cuffs
- 25 gsm sms fabric



Neck & Waist ties



Assurity Latex & Nitrile Gloves



- Premium quality
- Powder free
- Micro-textured
- Beaded cuff
- 200 per box Nitrile
- 100 per box Latex



XS S M L XL

Basic Dressing Pack

Order Code

TKBDP

- 1 x Sterile field
- 5 x Cotton balls
- 1 x Dressing tray
- 2 x Forceps Blue
- 1 x Dressing Towel
- 1 x Forceps Yellow
- 2 x Gauze swabs 7.5 x 7.5cm



Task Medical Tourniquet



25mm
Wide

Made in
Germany

- Premium quality
- Sturdy safety snap closure
- One-hand operation
- Latex free

Order Code
TKTOUR

Head Set Covers

Order Code

HEA

- Elastic rim
- Disposable
- Pack of 250
- Hygienic
- Sound and volume will not be affected



CT Biopsy Grid

Order Code

TKGRID

- Fast & accurate
- Easy to use
- Tape strips assure non-slipping positioning
- 1 cm markings



Finger Tip Pulse Oximeter

Order Code

TKPO300

- 4 direction display
- Large bright digital numbers
- Low power consumption
- User-friendly operation
- Automatically powers off



Faceshield

Order Code

FACESHMG

- Anti fog
- Splash resistant
- Light-weight
- Optically clear
- No dizziness



Multi Wipes



- 30 X 33cm
- Soft & strong
- 100% cellulose material
- Highly absorbent
- Dry wipes

Order Code
TKEW50

Disposable Prep Razor-SingleEdge

- Single Stainless Steel blade with comb
- Fixed type head
- Hard plastic handle
- Plastic cap to protect blade

Swedish Steel Blade
made from Europe's
most high quality steel.

Order Code
TKRAZ



Bite Block Covers



300 Pieces
per box

100 Boxes
per carton

- Latex free Polypropylene sleeves
- Easy application and removal
- Fully disposable

Order Code
TKBBC

Underpads

"blueys"

- Cost effective
- Highly absorbent
- High wicking
- Fibrous pulp
- Better than 5 ply equivalent
- Leak prevention



*Available in
2 weights*

Order Code

TKUP22

Order Code

TKUP30

Headpads

Order Code

TKHP50


- 3 Ply
- 50cm X 30cm
- Hygienic barrier
- Pillow protection
- Can also be used as a tray mat



Task Multi Roll

Available in 3 sizes

- | | | |
|-----------------------|----------------------|----------------------|
| ● Highly absorbent | ● Latex free | ● Polly backed |
| ● 2 ply perforated | ● Easy tear-off | ● Water resistant |
| ● 25cm x 50m | ● 50cm x 50m | ● 58cm x 50m |
| ● 12 rolls per carton | ● 6 rolls per carton | ● 6 rolls per carton |



25CM X 50M

TKMR25

50CM X 50M

TKMR50

58CM X 50M

TKMR58



Task Bed Sheet Roll

- 1 Ply
- 55cm x 80m
- 6 rolls per carton
- Added protection
- Keep furniture & surfaces clean between patients



Order Code

TKBR55

Pillow Sleeves

Blue

PILBLU

- Disposable
- 70cm x 45cm
- Water and blood resistant
- Open both ends



White

PILWHT

Pillow Protectors

Order Code

TKPILP

- Heavy duty vinyl
- 75cm x 50cm
- Easy to clean
- Zippered case
- Embossed



Fitted Bed Sheet



White

TKFITSWHT

- Non woven material
- Elastic edged fitted corners
- 200 x 75cm
- Disposable

Blue

TKFITSNB

Stretcher Sheet

White

SHEETWHT

- Strong, absorbent and water resistant
- 24 gsm paper
- 72cm width
- Cost effective

Blue

SHEET

Task Modesty Sheet

- XL Size 120cm x 80cm
- Soft 3 Ply paper
- Disposable
- Clean & hygienic
- Use as a small undersheet or to cover the patient



Order Code

TKMS80120

Vinyl Fitted Cover

- Environmental Touch Skin PU
- Water and oil resistant
- Elastic corners
- Fit both round and square table corners



Order Code

TKFITVN

Task Medical Premium

sms^{40gsm}
non woven fabric

Task Medical premium range consumables are made from a superior 40gsm SMS non woven fabric, comprising of two spunbond layers interlaid with one meltblown layer. They are an effective hygienic barrier against contamination on examination tables protecting your patients and furniture surfaces. Tear resistant and breathable.

Fitted bed Sheets



- 200cm x 70cm x 10m
- Disposable
- Fitted elastic ends
- Fitted corners
- Fluid resistant

Light Blue

TKPRFITS-LB

Modesty Sheets



- 100cm x 70cm
- Disposable
- Patient cover
- General purpose
- Fluid resistant

Light Blue

TKPRMODS-LB

Pillow cases



- 70cm x 45cm
- Disposable
- Sealed on 3 sides
- Open end flap
- Fluid resistant

Light Blue

TKPRPLCS-LB



CLARITY

Get to know us!
***We offer a wide range
of ultrasound products***

Used and preferred by sonographers
throughout Australia

Latex Probe Covers



**Non
Lubricated
TKPC**

**Lubricated
TKPCL**



Positioning Pad




- Fabric is breathable
- Antibacterial
- Waterproof
- Extreme durability
- Zipping water repellent
- Clean with standard clinical wipes

**Order Code
TKPOSTV**

ECG Electrode

**Order Code
CLEL3648**



- Fingerlift
- Sensor - Ag/AgCl
- Stud - Stainless steel
- Conductive gel - solid gel
- Carrier - PE foam 1mm thick
- Cover one side siliconized PET cover

Freezer Bags



**Order Code
CLFBW**

- Latex Free
- Single use only
- Individually Packed
- Size 45cm long x 30cm Wide
- Two Latex Free Elastic Bands
- Translucent Plastic Freezer Bags

**Includes 2
Elastic Bands**



Probe Cleaning Solution

- Prevents the transmission of bacteria from patient to patient
- It does not contain alcohol and aldehydes
- It has a handy spray cap
- It clears the image by removing gel residue on the probe

Odourless

**Order Code
CLPCS**

250ML



Ultrasound Transmission Gel

Clarity Ultrasound Gel ensures a clear and uninterrupted vision of the ultrasound waves on the screen by preventing the air flow between the skin and probe during all ultrasonography and doppler operations.

- Hypoallergenic, non-irritating
- Water soluble
- Non-staining easily cleanable
- Free from formaldehyde and salt
- No toxic effects
- Produced from harmless material
- Odourless
- Vacuum treated production
- Safe to use with the probe

**250ML
CLUSG250**

250ML



**5 Litre
CLUSG5L**



5 Litre

Sterile Ultrasound Gel

- It is used for ultrasonography and biopsy sampling processes in vitro fertilization, surgeries and vaginal and intra uterine examinations.
- Sterile pouch

**Order Code
CLSTUSG**

Hypoallergenic

20ML



Pack of 48

We've got you covered!

All of your dental needs in one place

TASK
MEDICAL



Dental Bib 2+1 & 1+1

All colours and patterns are available in 2+1 and 1+1



Saliva Ejectors



Head Rest Cover



Dental Roll 30 x 50m



Dental Cups 180ml

Now available in a range of colours




Dental Cotton Balls 300gr





Disposable & Reusable Holloware



Medicine Measures 30ml & 60ml



**Medicine
Measure Lid**

Clear

MML30

30ML

Clear

MM30

60ML

Clear

MM60

Gallipots 60ml & 150ml



60ML

Clear

BG60C

150ML

Clear

BG150C

Denture Bowls & Lids 250ml

Clear

B250CL50

Blue

B250BL50

Green

B250GL50



Kidney Dishes 700ml & 1000ml

Disposable & Autoclavable Available



Blue

Green

Clear

Yellow

Disposable Procedure Tray 200ml, 500ml & 1500ml



Blue

Green

Clear

Yellow

Disposable Bowls

Available in all sizes in color Clear, Blue & Green

250ml

500ml

1000ml

2000ml

5000ml

6000ml



Reusable Bowls

Available in all sizes in color Green & Yellow

150ml

500ml

1000ml

1500ml

2500ml

3500ml

6000ml



From nature to nurture

STAR ECO

Ecofriendly products



*Recyclable,
Biodegradable
& Compostable
products which
are disposable
and ecofriendly.*



Multi Cups 30ML



30ML
SPEC30

Recyclable
Latex Free
Cost effective
Liquid & Food Safe
Single Use
Light Weight
Multi-Functional



Multi Cups 50ML



Graduated Measuring Jugs



1000ML
SEMJ1000

General Purpose Semi-Bowls



Box of 100 **1500ML**
SEGPB1500

General Purpose Bowls



Box of 100

Male Urinals



1000ML
SEMU1000

Slipperpan Liners








Box of 100 **1300ML**
SESPL1300


Bedpan Liners



Box of 100






Injection Trays

-  Water and fluid resistant for up to 12 hours
-  Biodegradable and compostable within 12 weeks
-  500 Pieces/carton
-  200 mL capacity
-  Oil resistant



Box of 500 **200ML**
SEST200

Kidney Dishes

-  Water and fluid resistant for up to 12 hours
-  Biodegradable and compostable within 12 weeks
-  300 Pieces/carton
-  700 mL capacity
-  Oil resistant



Clinical Furniture

TASK
MEDICAL



Relieve back stress



- Improved posture
- Improved comfort
- Better ergonomics
- Premium Brand
- Lower Prices

WINBEX

World Class Quality

Authorized Australian Distributor



Mica Saddle Seat

Tilting seat, easy to change height position by using handle

Navy Blue
WX5113201NB

Black
WX113201BL

Premium Round Stool



Navy Blue
WX002327NB

Black
WX002327BL

Grey Teal
WX002327GT

Single lever for height adjustment

Comfort Line Stool




Navy Blue
WX002908NB

Black
WX002908BL

Grey Teal
WX002908GT

140mm / 200mm lift adjustment

Ascot Standard Saddle Seat



Design and built for ergonomics and comfort

Navy Blue
WX5105S02NB

Black
W51505S02BL

Kimi Saddle Seat



Helps to maintain correct posture

Navy Blue
WX511201NB

Black
WX511201BL

Doctors Stool T Seat



Comfort seat and backrest supports body and spine

Black
WX510710BL

Duma Doctor's Seat



Fully supported body and spine in a sedentary position

Navy Blue
WXDUMANB

Ascot Standard Saddle Seat Back



Backrest provides more support for longer working times

Navy Blue
WX511114NB

Black
WX511114BL

Kimi Saddle Seat With Back




Relieve stress from knees and hips

Navy Blue
WX511202NB

Black
WX511202BL

Microscope Chair



Separate adjustment for seat height, seat tilt and back angle

Black
WX5113301BL

Stainless Steel Trolleys

Task Medical instrument trolleys are manufactured using Hospital Grade 304 Stainless Steel. The trolley is durable whilst providing an aesthetic finish.




Side rails on all three sides, top and bottom



95mm Smooth running castors

Long Lasting & Strong

- High grade 304 stainless steel throughout
- Soft closing drawers
- German-made ball-bearing drawer runners
- 125mm deep drawers
- Locks on front castors



Powder Coated Trolley


Task Medical Powder Coated 1 & 2 Drawer Trolleys are classified as a high grade trolley with smooth castors and an outstanding finish. It is a perfect alternative to the stainless steel models and comes with a bottom shelf.

Durable & Strong

- Carbon steel powder coated
- Soft closing drawers
- Outstanding finish
- White colour



95mm Smooth running castors




Emergency Trolley


Task Medical Emergency Trolley has been specifically designed for the healthcare industry. This updated model is sturdy and with its larger wheels it allows for better mobility, highly responsive and effortless, when moved around.

- ABS moulded table with 304 stainless steel plate
- Four-sided raised stainless steel rim to avoid any falling of equipment
- High-quality aluminium moulded pillar
- Sideboard and backboard (Designed to be detachable)
- All drawers are detachable design and easy to clean
- Smart accessories
- Defib shelf
- Oxygen tank holder
- Central key lock
- Wire basket
- Extension table

Order Code
TKEMTAL5R



This trolley offers everything you need in an emergency at your fingertips - fast!



Treatment Trolley

Task Medical Treatment Trolley has a ABS moulded table with 304 stainless steel plate. The trolley has four-sided raised stainless rim's to help avoid any equipment falling.

- Aluminium drawers, side and back panels with aluminium pillars
- Soft closing removable drawers with inner dividers & ABS push handle
- Extended pull out side table
- 304 Stainless steel side basket
- Stainless Steel bench and shelf tops

Noise-free Castors With Front Brakes



Order Code
TK1DTRTWHT



Order Code
TK2DTRTWHT



Available with 1 or 2 Drawers in White colour

Blood Collection Chair

Task Medical Blood Collection Chair has been custom manufactured and designed to be the most comfortable and supportive chair.

- Premium finish that is durable, easy to use & clean
- Can also be used for examination
- Parents can sit with little ones & comfort during procedures
- Comfortable & supportive with a high-density soft foam cushion
- Integrated step stool
- Extra wide seat: 55cm

Black
TKBLDCHBL



Seated Capacity 200kg

Foldable Footrest



Grey
TKBLDCHGY

Mayo Instrument Trolley

The Task Medical Mayo Trolley is made from stainless steel construction. Some of the key components are its removable tray, adjustable height and lockable castors.

Height Adjustable

- Brush finished 304 stainless steel
- The top tray is detachable
- Four castors with brake
- Tray size 400 x 600mm
- Carrying capacity 35kg
- Dia 75mm castors

Order Code
TKMAYOTR



Linen Skip

Task Medical single linen skip trolleys are made from hospital grade tubular stainless steel. The linen skip trolley has four easy-glide casters with the front two castors featuring toe-snap brake locks.

Hospital Grade Stainless Steel

- Tubular stainless steel frame
- Front toe-snap castor locks
- Elevated bag support platform
- Size 38 x 54 x 85cm

75mm
Smooth
running castors



Order Code
TKLIN



Privacy Screen

This mobile Hospital Screen provides privacy for medical staff to examine patients. When not required, the Privacy Screen can be easily folded and stored.

Portable Privacy Screen

- Easy to clean
- Easily stored
- Open dimension: W180 x H185cm
- Panel width: 30cm/per panel
- Wheels for easy movement

Aluminium
Molded

Order Code
TKHSP



Step Stool

Task Medical Step Stools are durable and sturdy whilst being useful and efficient when extra height is required.

Heavy Duty

- Chromed steel
- Manufactured with 1" steel tubing
- Size 36 x 28 x 22cm
- 150kg Capacity

SILVER
VEIN
COLOR

Order Code
TKEST01F



Overbed Table

The Task Medical Overbed Table is practical and useful, by supporting patients with limited mobility to perform tasks such as eating, drinking, writing and much more.

Gas Strut Adjustable Height

- Top is made of ABS/PP
- Round corner design
- Powder coated iron base with plastic cover
- Silent casters with 2 brakes

Order Code
TKOBTABLE



Multi Purpose Mobility Bed

Task Medical EzyM Multipurpose Mobility Bed

A manual wind hi-lo patient mobility bed is ideal for any medical practice

The new EzyM Mobility Bed from Task Medical is a step up from the traditional heavy steel-framed, flatbed patient gurney. Almost entirely made from a high impact ABS polymer and a high strength alloy frame, this new generation of the gurney is both a modern patient transport stretcher and a hi-lo bed. With a multitude of uses that would be beneficial for service in any kind of medical practice, aged care facility or hospital, the EzyM would prove to be a sound investment given its low price.

EZYM HI LO
by TASK MEDICAL



Ideal for any medical practice!

Order Code
TKEZYM

- Recovery room
- Procedure room
- Patient transporter
- Infusion therapy bed
- Emergency stretcher
- 180kg SWL
- Manual wind
- 56-86cm



Includes new extra thick mattress upgrade

Hospital grade easy clean vinyl with 60mm thick comfort padding

two oxygen cylinder holder rings



Easy roller mri safe wheel chair

MR SAFE Certified MR safe (7 Tesla) by Siemens and Phillips

Made from recyclable PE and PA thermoplastics

100% non metallic non magnetic

EasyRoller has been especially designed as a MR SAFE patient transport device for use in hospitals and MR laboratories. Patients can easily be moved through the entire hospital and all the way through the MRI room without interfering with the MRI equipment. Norwegian design, German made.



Transport patients directly to the mri machine

MRI Safe Wheelchair



- No metal parts - Non metallic
- No maintenance recyclable plastic
- Easy maneuvering swivel wheel axle
- Tubeless soft ride tyres & park brake
- Double POM glass ball bearing wheel axles
- Waterproof padded seat, back & armrests
- Approved Standard EN12183 and ISO 7176

Order Code
WHEEASYRL




Camino Medical Hi-Lo Electric Exam Couch

With/without 360° walk-free foot bar control
3 section - Single Motor


This is a Medical Exam Couch with exceptional strength and stability which outperforms many in its class

- Foot pedal control
- Face hole
- Thicker 64mm foam padding
- Gas strut tilt control
- Steel frame
- Front bar key lock
- Lower & raise retractable castors
- Adjustable feet
- Fireproof PU/leather
- 250kg SWL

Black
102531
Charcoal
102531
Agate Blue
102360



Walk free control bar Agate Blue
102361



Superior German made motor

WARRANTY
5 YEAR FRAME
3 YEAR ELECTRICAL



250kg SWL

Camino Physio Hi-Lo Electric Treatment Table

Designed to support multiple posture adjustments for physio treatments

- Electric height adjust.
- 46-95cm. SWL 250kg
- A double-sided synchronous pneumatic control of the tail-adjustment mechanism makes it more sturdy, and achieves quiet and fast adjustment

Dual struts



Lower & raise retractable castors



Longer section


Agate Blue
102535

Sonora All Electric Medical Aesthetic Chair

Beauty & Cosmetic
3 section - 3 Motor
Available in Black or White

The Sonora contributes to serve beauty customers with more comfort and indulgence and to improve the treatment capabilities and convenience of beauticians and medical aesthetic clinicians.

- Electric height adjust.
- 56-90cm. SWL 200kg
- Electric backrest 0°-75° angle adjustment
- Electric leg rest 0°-90° angle adjustment
- Body Oil resistant PU leather
- Luxurious 82mm thick high density
- CA-117 fireproof foam padding
- Removable leveled side armrests
- Removable /adjust' neck/face pillow



3 Motors
Height motor
Backrest motor
Footrest motor

Face hole

180° swing back arms

Superior German made Motors

More positions

Height & multi-position hand controller

Stylish mono frame lift mechanism

Sturdy base

3 Motor control

Adjustable foot rest extension

Black 102363
White 102362



Standard Examination Couch

**First Class quality desing & construction
affordably priced!**

This is a heavy duty examination couch made from heavy grade steel providing an exceptional sturdy base an platform for your patients.

- Australian Desing
- Easy hand lift backrest with 10 stage angle adjustment hinge
- Heavy Duty round tubular legs
- Heavy grade steel constructed frame
- Foam Density: 36g/m2
- Heavy Duty PVC vinyl
- Easy clean and hospital grade
- 260kg SWL

**Available
Black, Grey & Navy Blue**



Standard Fixed Height Examination Couch

- Fixed height 82cm
- Two section
- Adjustable angle, lift up back rest
- Powder coated steel frame
- New easy lift and lock mechanism for easy 2 minute assembly
- Comes flat packed, easy to assemble
- Padding thickness 5cm, Foam Density: 36g/m2
- Total Length: 195cm
- Width: 71cm
- Height: 82cm
- 250kg SWL

250kg SWL



MRI Safe Stretcher

With adjustable height and drip stand

The Task Medical MRI Conditional stretcher is Non Magnetic and compatible / suitable for 1.5T and 3.0T MR conditional equipment.

- Aluminium construction
- Adjustable (hand crack) height*
- Adjustable back rest
- Safe Working Load 169kg
- 70mm thick padded mattress
- Drop down side rails
- Stainless steel IV pole
- 125mm castors, 2 with brakes
- Back rest adjustment range 0-75 degrees (lever)
- Max Height 189cm
- Min Height 46cm
- Overall Length 189cm
- Overall Width 56cm



Multipurpose Treatment Chair

Hi-Lo 3 Section / 2 Motor / Grey Upholstery

Task Medical's podiatry chair focuses on patient comfort, with armrests and a long backrest reaching head height. The position variations provide better positions for the practitioner's examination or treatment procedures.

- Double gas strut back adjustment and single gas strut foot adjustment,
- seat and height of table adjustable by actuators
- Four retractable medical grade castors
- Four large pad levelling feet
- 1 pair adjustable tilt back armrest
- Medical grade PVC vinyl
- Size 190mm L x 600mm W
- Adjustable height range 510-950mm
- Adjustable backrest 0 degrees to +78 degrees
- Adjustable footrest 0 degrees to -90 degrees
- 250kg SWL

**250kg
SWL**



Deluxe Gynaecology Chair With Stirrup & Foot Boards

3 Section / 2 Motor / Deluxe Gynaecology Chair - Table / Blue

- Double gas strut back adjustment
- Four retractable medical grade castors
- 2 large rear pad levelling feet
- 1 pair adjustable arm rests
- 1 pair stirrup leg supports
- 1 pair foot boards
- 1 stainless steel debris tray
- Medical grade PVC vinyl upholstery
- Size 1820 x 600mm
- Adjustable back: 0 to 80 degrees
- Adjustable seat: 0 to +20 degrees
- Adjustable height: 590-1060mm
- 2 piece actuator with foot switch adjustable height and seat\
- 250kg SWL

**250kg
SWL**



Stainless Steel Mobile IV Stand With 4 Hooks & Chrome Base

Mobile IV Stand with 4 hooks

- Mobile support pole for intravenous use
- Adjustable height
- Double hooks for dual capacity
- Five star metal base for maximum stability
- Swivel castors
- Durable chromed steel pole

TASK
MEDICAL

