


STAR ECO

from nature to nurture

Star Eco introduces its comprehensive range of disposable and ecofriendly products.



1000mL Graduated Measuring Jugs

REF SEMJ1000

100 pieces per carton



- Biodegradable and compostable within 12 weeks (test report available) including carton packaging
- Maceratorable and recyclable*
- Made with agro fibre additive technology offering significant environmental and sustainability advantages and lifestyle improvements to farmers through income supplementation
- Directs agricultural bi-products into important raw materials
- Water and fluid stress resistant for up to 12 hours. Products comply fully with the internationally recognized PAS 29:1999 standards (ability to hold water/liquid at 38°C for a minimum of 4 hours)

- Oil resistant
- Non toxic – free from heavy metals and pesticides and processed via natural means without bleach or other chemicals
- Latex and wax free
- Non sterile and Single Patient Use for Contamination Control
- Puncture resistant
- Ordourless
- Low fibre shred
- Rigid, smooth and durable finish
- Stackable to reduce storage space requirements
- Offers storage stability advantages
- Cost effective and economical

* Recyclable after correct decontamination procedures and guidelines.

	Materials	Appearance	Typical Product Thickness
Star Eco	Recycled corrugated cartons mixed with agro fibre materials	Rigid, smooth and durable	0.8 to 1.0 mm
Others	Old, unclean corrugated cartons mixed with office paper and bagasse	Uneven, flimsy and easily folded	0.4 to 0.6 mm

Ideal for use in General Practices, Dental Practices, Medical Clinics, Aged Care Facilities and Hospitals.