

STAR ECO

from nature to nurture

Star Eco introduces its comprehensive range of disposable and ecofriendly products.



1500mL General Purpose Semi-Bowls

REF SEGPB1500

100 pieces per carton



- 🌿 Biodegradable and compostable within 12 weeks (test report available) including carton packaging
- 🌿 Maceratorable and recyclable*
- 🌿 Made with agro fibre additive technology offering significant environmental and sustainability advantages and lifestyle improvements to farmers through income supplementation
- 🌿 Directs agricultural bi-products into important raw materials
- 🌿 Water and fluid stress resistant for up to 12 hours. Products comply fully with the internationally recognized PAS 29:1999 standards (ability to hold water/liquid at 38°C for a minimum of 4 hours)

- 🌿 Oil resistant
- 🌿 Non toxic – free from heavy metals and pesticides and processed via natural means without bleach or other chemicals
- 🌿 Latex and wax free
- 🌿 Non sterile and Single Patient Use for Contamination Control
- 🌿 Puncture resistant
- 🌿 Odourless
- 🌿 Low fibre shred
- 🌿 Rigid, smooth and durable finish
- 🌿 Stackable to reduce storage space requirements
- 🌿 Offers storage stability advantages
- 🌿 Cost effective and economical

* Recyclable after correct decontamination procedures and guidelines.

| | Materials | Appearance | Typical Product Thickness |
|-----------------|---|----------------------------------|---------------------------|
| Star Eco | Recycled corrugated cartons mixed with agro fibre materials | Rigid, smooth and durable | 0.8 to 1.0 mm |
| Others | Old, unclean corrugated cartons mixed with office paper and bagasse | Uneven, flimsy and easily folded | 0.4 to 0.6 mm |

Ideal for use in General Practices, Dental Practices, Medical Clinics, Aged Care Facilities and Hospitals.