

# STAR ECO

from nature to nurture

Star Eco introduces its comprehensive range of disposable and ecofriendly products.



## 3000mL General Purpose Bowls

REF SEGPB3000

100 pieces per carton



- 🌿 Biodegradable and compostable within 12 weeks (test report available) including carton packaging
- 🌿 Maceratorable and recyclable\*
- 🌿 Made with agro fibre additive technology offering significant environmental and sustainability advantages and lifestyle improvements to farmers through income supplementation
- 🌿 Directs agricultural bi-products into important raw materials
- 🌿 Water and fluid stress resistant for up to 12 hours. Products comply fully with the internationally recognized PAS 29:1999 standards (ability to hold water/liquid at 38°C for a minimum of 4 hours)

- 🌿 Oil resistant
- 🌿 Non toxic – free from heavy metals and pesticides and processed via natural means without bleach or other chemicals
- 🌿 Latex and wax free
- 🌿 Non sterile and Single Patient Use for Contamination Control
- 🌿 Puncture resistant
- 🌿 Odourless
- 🌿 Low fibre shred
- 🌿 Rigid, smooth and durable finish
- 🌿 Stackable to reduce storage space requirements
- 🌿 Offers storage stability advantages
- 🌿 Cost effective and economical

\* Recyclable after correct decontamination procedures and guidelines.

	Materials	Appearance	Typical Product Thickness
<b>Star Eco</b>	Recycled corrugated cartons mixed with agro fibre materials	Rigid, smooth and durable	0.8 to 1.0 mm
<b>Others</b>	Old, unclean corrugated cartons mixed with office paper and bagasse	Uneven, flimsy and easily folded	0.4 to 0.6 mm

**Ideal for use in General Practices, Dental Practices, Medical Clinics, Aged Care Facilities and Hospitals.**

Distributed by:  
Constar X Pty Ltd



Unit 4, 68 Parramatta Road, Underwood QLD 4119

📞 1300 778 924

✉ info@constar.com.au



constar.com.au/star-eco