



CLARITY

Ultrasound Gels, Probe Covers & Probe Cleaning Solutions for performing ultrasonography procedures.



Clarity Ultrasound Gel 250ml and 5L

- Not affected by salt on the skin during usage
- Ensures a comfortable ultrasonography procedure
- Hypoallergenic, does not irritate the skin
- Doesn't contain oil or oily agents, formaldehyde or salt
- No toxic effect or odour



Clarity Sterile Ultrasound Gel 20ml

- Sterilised with gamma rays, this gel is used for ultrasonography and biopsy sampling procedures, in vitro fertilization, surgeries and vaginal and intra uterine examinations.



CLSTUSG

Clarity Probe Covers for Transvaginal Ultrasound

- Specially designed for Gynaecological Ultrasonography. Each probe cover is rolled for easy application to cover a variety of trans-vaginal transducers. 100 individually foil sealed probe covers in each dispenser box.
- Designed and manufactured specifically to meet transducer probe cover specifications
- Manufactured from high quality homogenised low allergy rubber latex
- Meets Quality Standard ISO 4070
- Length: 200mm Circumference: 53mm Width: 34mm
- Lubricated and Non-Lubricated probes available
- Single use / Non Sterile



Clarity Probe Cleaning Solution 250ml

- Ultraseptic ultrasound probe cleaning spray solution removes residual gel left over the probe.
- Helps prevent contamination from patient to patient thanks to its active agents.



CLPCS