


equipment &
consumables

TASK
MEDICAL



taskmedical.com.au



X ray gown

- Belt loop
- Disposable
- Shoulder sleeve
- Generous size
- Box of 100

Standard **GOWNT** X large **GOWNXL**

Isolation gown

- Belt loop
- Long sleeve
- Disposable
- Box of 100

Standard **GOWNY**



Face mask

- High filtration
- 3 ply
- Ear loops
- Pack of 50



TKO301

Task Medical Tourniquet

Quality product made in Germany

Proven for daily use in hospitals, medical practices, and for paramedics.

- Skin-friendly, flexible and soft ribbon 25 mm wide
- Sturdy safety snap closure
- One-hand operation
- Metalfree
- No squeezing or pinching of the skin
- Latex free

TKTOUR



Head set covers

- Elastic rim
- Disposable
- Pack of 250
- Hygienic

HEA



Fitted couch cover

- Non woven
- Disposable
- 200cm x 75cm
- Carton of 100

White

TKFITWHT

Navy Blue

TKFITSNB



Stretcher sheet

- Non woven
- 240cm x 70cm
- Disposable
- Carton of 100

White

SHEETWHT

Navy Blue

SHEET



Pillow sleeves

- Disposable
- 70cm x 45cm
- Disposable
- Carton of 200

White

PILWHT

Navy Blue

PILBLU



TASK MEDICAL PREMIUM

sms 40gsm
non woven fabric

Task Medical premium range consumables are made from a superior 40gsm SMS non woven fabric, comprising of two spunbond layers interlaid with one meltblown layer. They are an effective hygienic barrier against contamination on examination tables protecting your patients and furniture surfaces. Tear resistant and breathable.

Fitted bed sheets

- 200cm x 70cm x 10 cm
- Disposable
- Fitted elastic ends
- Fitted corners
- Fluid resistant

Light blue

TKPRFITS-LB



Modesty sheets

- 100cm x 70cm
- Disposable
- Patient cover
- General purpose
- Fluid resistant

Light blue

TKPRMODS-LB



Pillow cases

- 70cm x 45 cm
- Disposable
- Sealed on 3 sides
- Open end flap
- Fluid resistant

Light blue

TKPRPLCS-LB



Versatile, popular product



TASK Multi Roll Available in 3 sizes

- Highly absorbent
- 2ply Perforated
- 25cm x 50m
- 12 rolls per carton

- Highly absorbent
- 2ply Perforated
- 50cm x 50m
- 6 rolls per carton

- Highly absorbent
- 2ply Perforated
- 58cm x 50m
- 6 rolls per carton



TKMR25

TKMR50

TKMR58

TASK Bed Roll


- Keep furniture & surfaces clean between patients
- Added protection
- 1 ply
- 55cm x 80m
- 6 rolls per carton



TKBR55

TASK Utility Roll - Laminated

- Fluid proof & hygienic barrier
- Perforated Polyethylene Backed
- 68cm x 50m
- 12 rolls per carton



TKUR6850

TASK Modesty sheet

- 120cm x 80cm
- 3ply Disposable
- Case of 100
- Protect & Respect



XL Size
TKMS80120

Basic dressing pack

- 1 x Sterile field
- 1 x Dressing tray
- 1 x Dressing towel
- 5 x Cotton balls
- 2 x Forceps blue
- 1 x Forceps yellow
- 2 x Gauze swabs 7.5cm x 7.5cm



TKBDP **QUALITY** Tested *****
STERILE EO 0197

Underpads


TASK underpads are ideal for fluid management, providing extra protection for beds, trolleys, trays etc.

- High wicking
- Fibrous pulp
- Leak prevention
- Highly absorbent
- Sealed 4 sides

Premium
60cm x 40cm
Ctn 300
TKUP30

Economy
60cm x 40cm
Ctn 300
TKUP22

Half Size
Economy
30cm x 40cm
Ctn 400
5 ply
TKUP25



Everyday wipes

- Highly absorbent
- Dry wipes
- Latex free
- 30cm x 33cm
- 50 wipes / pack
- Carton of 12 packs



TKEW50

Head pads

- Hygienic barrier
- Pillow protection
- 50cm x 30cm
- 3 ply
- Carton of 1000



TKHP50

Relieve back stress



- ▣ Improved posture
- ▣ Improved comfort
- ▣ Better ergonomics
- ▣ Premium Brand
- ▣ Lower Prices



Mica Saddle Seat

Tilting seat, easy to change height position by using handle

WINBEX
World Class Quality

Authorised Australian Distributor

Navy Blue
WX5113201NB

Black
WX113201BL

Premium quality stools at affordable prices



Premium Round Stool

Single lever for height adjustment

Navy Blue
WX002327NB
Black
WX002327BL
Grey Teal
WX002327GT



Comfort Line Stool

140mm / 200mm lift adjustment

Navy Blue
WX002908NB
Black
WX002908BL
Grey Teal
WX002908GT



Ascot Standard Saddle Seat

Design and built for ergonomics and comfort

Navy Blue
WX5105S02NB
Black
WX5105S02BL



Kimi Saddle Seat

Helps to maintain correct posture

Navy Blue
WX511201NB
Black
WX511201BL



Doctors Stool T Seat

Comfort seat and backrest supports body and spine

Black
WX510710BL



Duma Doctor's Seat

Fully supported body and spine in a sedentary position

Navy Blue
WXDUMANB



Ascot Standard Saddle Seat Back

Backrest provides more support for longer working times


Navy Blue
WX511114NB
Black
WX511114BL



Kimi Saddle Seat with Back

Relieve stress from knees and hips

Navy Blue
WX511202NB
Black
WX511202BL



Microscope Chair

Separate adjustment for seat height, seat tilt and back angle


Black
WX5113301BL



e200 transport stretcher

Equipped with a four-wheel brake system and a central fifth wheel which offers excellent mobility to provide the best safety for patients, the e200 provides durable quality and a 300kg weight capacity.

TKST200



- Foot pedal for hydraulic system control on both sides
- Integrated storage tray
- Two-section bed board / Collapsible side rails
- Drainage Bag clasp x 2 / Urinary Bag clasp x 2
- Foot controlled brake / Retractable steering fifth wheel
- Steer pedal for four wheels at both ends
- Head-end delivery push bar / Foot-end assist push bar
- Dual pneumatic assisted backrest
- Two stage movable IV pole x 1
- Two stage permanent folding IV pole x 1
- Velcro tape underneath the mattress
- Protective corner bumpers x 4
- Hook / jack for upright oxygen tank holder
- 300 kg SWL

Multi Purpose Mobility Bed

EZYM HI LO
by TASK MEDICAL

- Recovery room
- Procedure room
- Patient transporter
- Infusion therapy bed
- Emergency stretcher
- 180kg SWL
- Manual wind
- 56-86cm

Ideal for any medical practice!

TKEZYM



Easy roller mri safe wheel chair

MR SAFE Certified MR safe (7 Tesla) by Siemens and Phillips

Made from recyclable PE and PA thermoplastics

100% non metallic non magnetic

EasyRoller has been especially designed as a MR SAFE patient transport device for use in hospitals and MR laboratories. Patients can easily be moved through the entire hospital and all the way through the MRI room without interfering with the MRI equipment. Norwegian design, German made.



- No metal parts - Non metallic
- No maintenance recyclable plastic
- Easy maneuvering swivel wheel axle
- Tubeless soft ride tyres & park brake
- Double POM glass ball bearing wheel axles
- Waterproof padded seat, back & armrests
- Approved Standard EN12183 and ISO 7176



Transport patients directly to the mri machine


MRI Safe Wheelchair

WHEEASYRL



TASK
MEDICAL

taskmedical.com.au



Single no drawer



TK0D50N

Single one drawer



TK1D50N

Single two drawer



TK2D50N

Single three drawer



TK3D50N

Single four drawer



TK4D50N

Mayo instrument



TKMT3L

Stainless steel trolleys

Task Medical trolleys are made from high grade 304 stainless steel throughout, with double safety folds and smooth welding.

The top and bottom shelves are surrounded by safety rails on all three sides.

Movability is made easy by 95mm smooth rolling castors with lock brakes on the front.

The drawers are soft closing magnetic with German made ball bearing runners.



TK0D80N



TK2D80N



TK4D80N

Emergency Trolley

Have your CPR & Resus' gear ready to go in the one unit

TKEM65032



Clinical Trolley

TKCT1D

Extra large castors for effortless maneuverability



Camino Medical Hi-Lo Electric Exam Couch

With/without 360° walk-free foot bar control

3 section - Single Motor




Superior German made motor

250kg SWL

WARRANTY
5 YEAR FRAME
3 YEAR ELECTRICAL

This is a Medical Exam Couch with exceptional strength and stability which outperforms many in its class

- Foot pedal control
- Face hole
- Thicker 64mm foam padding
- Gas strut tilt control
- Steel frame
- Front bar key lock
- Lower & raise retractable castors
- Adjustable feet
- Fireproof PU/leather
- 250kg SWL



Black
102531

Charcoal
102531

Agate blue
102360

Walk free control bar
Agate blue
102361



Camino Physio Hi-Lo Electric Treatment Table

Designed to support multiple posture adjustments for physio treatments

- Electric height adjust.
- 46-95cm. SWL 250kg
- A double-sided synchronous pneumatic control of the tail-adjustment mechanism makes it more sturdy, and achieves quiet and fast adjustment

Dual struts



Lower & raise retractable castors

Longer section

Agate blue
102535



Sonora All Electric Medical Aesthetic Chair


Beauty & Cosmetic

3 section - 3 Motor

Available in Black or White

The Sonora contributes to serve beauty customers with more comfort and indulgence and to improve the treatment capabilities and convenience of beauticians and medical aesthetic clinicians.

- Electric height adjust.
- 56-90cm. SWL 200kg
- Electric backrest 0°-75° angle adjustment
- Electric leg rest 0°-90° angle adjustment
- Body Oil resistant PU leather
- Luxurious 82mm thick high density
- CA-117 fireproof foam padding
- Removable leveled side armrests
- Removable /adjust' neck/face pillow



3 Motors
Height motor
Backrest motor
Footrest motor

Height & multi-position hand controller

Stylish mono frame lift mechanism

Sturdy base

Black
102363

White
102362

Face hole

180° swing back arms

Superior German made Motors

More positions

3 Motor control

Adjustable foot rest extension



Dental Bib 2+1 & 1+1

All colours and patterns are available in 2+1 and 1+1



Saliva Ejectors



Head Rest Cover



Dental Roll 30cm x 50m



Dental Cups 180ml

Now available in a range of colours



Dental Cotton Rolls 300gr

