



CLARITY

Ultrasound Gels, Probe Covers & Probe Cleaning Solutions for performing ultrasonography procedures.



Clarity Ultrasound Gel 250ml and 5L

- Not affected by salt on the skin during usage
- Ensures a comfortable ultrasonography procedure
- Hypoallergenic, does not irritate the skin
- Doesn't contain oil or oily agents, formaldehyde or salt
- No toxic effect or odour



Clarity Sterile Ultrasound Gel 20ml

- Sterilized with gamma rays, this gel is used for ultrasonography and biopsy sampling procedures, in vitro fertilization, surgeries and vaginal and intra uterine examinations.



Clarity Probe Covers for Transvaginal Ultrasound

- Specially designed for Gynecological Ultrasonography. Each probe cover is rolled for easy application to cover a variety of trans-vaginal transducers. 100 individually foil sealed probe covers in each dispenser box.
- Designed and manufactured specifically to meet transducer probe cover specifications
- Manufactured from high quality homogenized low allergy rubber latex
- Meets Quality Standard ISO 4070
- Length: 200mm Circumference: 53mm Width: 34mm
- Lubricated and Non-Lubricated probes available
- Single use / Non Sterile



Cleaning Solution 250ml

- Ultraseptic ultrasound probe cleaning spray solution removes residual gel left over the probe.
- Helps prevent contamination from patient to patient thanks to its active agents.



Clarity Freezer Bags

- Single Use
- 30cm x 45cm
- 2 x latex free bands
- Individually packed



Bed Rolls and Multi Rolls

Our Bed rolls and multi rolls can be used as surface protective covers for exam couches and other medical furniture, general purpose wipes for cleaning, hand towels, tray liners and to soak up most fluid spills. Made from a non woven paper like fabric the rolls come in a range of sizes, perforated and non perforated.

Bed sheet roll non perforated

- ✔ Single Ply **TKBR55**
- ✔ 55cm x 80m
- ✔ Highly absorbent
- ✔ Latex free



Bed sheet roll perforated

- ✔ 2 Ply **TKMR58**
- ✔ 58cm x 80m
- ✔ Highly absorbent
- ✔ Perforated every 50cm



Bed sheet roll perforated polyethylene backed

- ✔ 1 Ply **TKUR68**
- ✔ 68cm x 50m
- ✔ Poly backed
- ✔ Perforated every 38cm



General purpose multi roll perforated

- ✔ 1 Ply **TKMR50**
- ✔ 50cm x 50m
- ✔ Highly absorbent
- ✔ Perforated every 50cm



Headpads

A hygienic barrier to treatment surfaces, keeping the head and face pad area of pillows and exam couches free of oils and perspiration.

- ✔ Protection for treatment surfaces
- ✔ Hygienic barrier
- ✔ 50cm x 30cm
- ✔ 3 Ply



TKHP50



Modesty Sheets

The use of disposable sheets removes issues associated with collection, washing and cleanliness if reusable linen.

- ✔ White
- ✔ 120cm x 80cm
- ✔ Disposable
- ✔ Non woven fabric
- ✔ Case of 100



TKMS80120

Everyday Wipes

Everyday wipes are a high-strength, super soft, and highly absorbent cleaning wipe designed for healthcare providers. Can safely be used on patients skin, medical equipment and surfaces.

- ✔ 50 Dry wipes
- ✔ 30cm x 33cm
- ✔ White
- ✔ 100% Cellulose material
- ✔ Biodegradable
- ✔ Latex Free
- ✔ Non-Sterile
- ✔ Single Use



TKEW50