

# assurity

ambidextrous  
textured finger tips

powder free

premium  
quality

## latex examination gloves

non sterile  
chlorinated

*when personal protection and patient care  
is in your hands - you want assurity*



# assurity

made for medical for the protection of healthcare professionals



extra small  
aslegnxs

small  
aslegns

medium  
aslegnm

large  
aslegnl

extra large  
aslegnxl

# assurity

## latex examination gloves



*glove conforms to*

- EN455 parts 1,2,3 & 4
- ASTM D3578
- ISO9001:2015
- ISO13485:2016



Assurity Nitrile Gloves can be purchased by the box or carton. There are ten boxes in a carton and 100 gloves in a box. So per carton, that's a total of 1000 gloves at a much lesser price.




bx100



*natural colour*


Assurity latex examination gloves offer healthcare professionals improved biological barrier protection. Powder free, the gloves are made from natural rubber latex. The robust middleweight glove offers consistent sensitivity and tensile strength. Assurity latex examination gloves provide long-lasting comfort, durability and decreased hand fatigue and improved dexterity and sensitivity. They conform and are tested to the highest international standards. Assurity Latex Examination Gloves are available in five different sizes with individually colour coded dispenser boxes.



- premium quality
- 4.5mil latex
- powder free
- natural rubber latex
- micro-textured
- ambidextrous
- non sterile single use
- beaded cuff
- AQL 1.5

ARTG ID 392396

### sizes & dimensions



|    |      |         |             |                |
|----|------|---------|-------------|----------------|
| xs | 70   | ± 10 mm | min. 240 mm | ① min. 0.13 mm |
| s  | 80   | ± 10 mm | min. 240 mm | ② min. 0.10 mm |
| m  | 95   | ± 10 mm | min. 240 mm | ③ min. 0.07 mm |
| l  | 110  | ± 10 mm | min. 240 mm |                |
| xl | ≥110 |         | min. 240 mm |                |

### single wall thickness

### physical properties

| assurity          | unaged        | aged          |
|-------------------|---------------|---------------|
| tensile strength  | min. 18.0 MPa | min. 14.0 MPa |
| elongation        | min. 650%     | min. 500%     |
| strength at break | min. 6N       | min. 6N       |

distributed by constar x pty ltd  
unit 4, 68 parramatta road, underwood qld 4119



1300 788 924  
07 5450 7092

info@constar.com.au  
www.constar.com.au