

assurity

ambidextrous
textured finger tips

powder free

premium
quality



latex examination gloves

non sterile
chlorinated

*when personal protection and patient care
is in your hands - you want assurity*



assurity

made for medical for the protection of healthcare professionals



extra small
aslegnxs

small
aslegns

medium
aslegnm

large
aslegnl

extra large
aslegnxl

constar.com.au



assurity

latex examination gloves



Assurity Nitrile Gloves can be purchased by the box or carton. There are ten boxes in a carton and 100 gloves in a box. So per carton, that's a total of 1000 gloves at a much lesser price.



bx100




single use

glove conforms to

- EN455 parts 1,2,3 & 4
- ASTM D3578
- ISO9001:2015
- ISO13485:2016


Assurity latex examination gloves offer healthcare professionals improved biological barrier protection. Powder free, the gloves are made from natural rubber latex. The robust middleweight glove offers consistent sensitivity and tensile strength. Assurity latex examination gloves provide long-lasting comfort, durability and decreased hand fatigue and improved dexterity and sensitivity. They conform and are tested to the highest international standards. Assurity Latex Examination Gloves are available in five different sizes with individually colour coded dispenser boxes.

- premium quality
- 4.5mil latex
- powder free
- natural rubber latex
- micro-textured
- ambidextrous
- non sterile single use
- beaded cuff
- AQL 1.5



ARTG ID 39236

sizes & dimensions



Size	Width (mm)	Length (mm)
xs	70 ± 10 mm	min. 240 mm
s	80 ± 10 mm	min. 240 mm
m	95 ± 10 mm	min. 240 mm
l	110 ± 10 mm	min. 240 mm
xl	≥110	min. 240 mm

single wall thickness



- ① min. 0.13 mm
- ② min. 0.10 mm
- ③ min. 0.07 mm

physical properties

assurity	unaged	aged
tensile strength	min. 18.0 MPa	min. 14.0 MPa
elongation	min. 650%	min. 500%
strength at break	min. 6N	min. 6N

distributed by constar x pty ltd
unit 4, 68 parramatta road, underwood qld 4119



1300 788 924
07 5450 7092

info@constar.com.au
www.constar.com.au