

assurity

latex free
powder free

ambidextrous
textured finger tips



premium
quality

non sterile
lab chemical tested

nitrile examination gloves

*when personal protection and patient care
is in your hands - you want assurity*

assurity



made for medical for the protection of healthcare professionals



extra small
asnegbxs

small
asnegbs

medium
asnegbm

large
asnegbl

extra large
asnegbxl

constar.com.au



assurity

nitrile examination gloves



glove conforms to

- EN455 parts 1,2,3 & 4
- ASTM D6319 & ASTM D6978
- ISO9001:2015
- ISO13485:2016




Assurity Nitrile Gloves can be purchased by the box or carton. There are ten boxes in a carton and 200 gloves in a box. So per carton, that's a total of 2000 gloves at a much lesser price.



blue


Assurity nitrile examination gloves offer improved biological barrier protection while mimicking the slim fit and comfort of a latex glove. Powder free, the gloves are made from 100% Nitrile (Non Latex). The lightest of our middleweight glove offers consistent sensitivity and tensile strength and its Chemotherapy and Fentanyl rating makes this a popular choice in healthcare facilities. Assurity nitrile examination gloves provide long-lasting comfort, durability and decreased hand fatigue. Assurity gloves conform and are tested to the highest international standards.



- premium quality
- 3.2mil powder free nitrile
- latex free
- textured finger tips
- ambidextrous
- non sterile single use
- beaded cuff
- chemotherapy drugs tested
- AQL 1.5


ARTG ID 39234

sizes & dimensions



xs	70 ± 10 mm	min. 240 mm
s	80 ± 10 mm	min. 240 mm
m	95 ± 10 mm	min. 240 mm
l	110 ± 10 mm	min. 240 mm
xl	≥110	min. 240 mm

single wall thickness



- ① min. 0.10 mm
- ② min. 0.06 mm
- ③ min. 0.05 mm

physical properties

assurity	unaged	aged
tensile strength	min. 14.0 MPa	min. 14.0 MPa
elongation	min. 500%	min. 400%
strength at break	min. 6N	min. 6N

*

distributed by constar x pty ltd
unit 4, 68 parramatta road, underwood qld 4119



1300 788 924
07 5450 7092

info@constar.com.au
www.constar.com.au