

Stainless Steel Trolley 4 Drawers

TASK
MEDICAL

Task Medical instrument trolleys are manufactured using Hospital Grade 304 Stainless Steel. The trolley is durable whilst providing an aesthetic finish.

Long Lasting & Strong

- High grade 304 stainless steel throughout
- Soft closing drawers
- German-made ball-bearing drawer runners
- Locks on front castors
- Smooth welding finish

Side rails on all
three sides



4 x 125mm
deep shelves



95mm
Smooth
running castors



www.taskmedical.com.au



SKU
TK4D50N