

# Fenestrated Drape 50cm x 68cm




Task Medical Fenestrated Drape, is ideal for many surgical and therapeutic procedures and to help protect patients and medical staff from bacterial transfer.

## Product Information

- Peel pack
- Antistatic
- 4 ply poly-lined
- Individually sterile packed
- Folded design that is easy to open
- 10cm diameter fenestration
- Non-woven fabric
- Single use only
- EO sterilised

50cm



68cm



## Packaging

- 20 pouches / box
- 18 boxes / carton
- Total pieces / carton - 360

### Carton Size

### Box Size

- 54x32x45cm
- G.W 11kg
- N.W 9kg
- 20x15x17cm
- Weight 0.50kg

CODE  
TKDRF5068



[www.constar.com.au](http://www.constar.com.au)

