

## Mixing Cannula 14cm




Task Medical Mixing Cannula, is a semi-rigid sterile tube designed to facilitate accurate and contamination-free mixing of liquids or solutions in healthcare, laboratory and research environments.

### Product Information

- Peel pack
- Precision reach and manoeuvrability
- Compatible with various solutions or fluids
- User-friendly design with a tapered tip
- Secure attachment mechanism
- Ideal for sterile mixing and chemical applications
- Complies with medical industry standards
- Single-use for infection control
- Designed to be attached to standard luer slip syringes
- EO sterilised

*14 cm long*



### Packaging

- 100 pieces / bag

### Carton Size

- 41x28x18.5cm
- G.W 3.05kg
- N.W 2.55kg

CODE  
TKMIX



[www.constar.com.au](http://www.constar.com.au)

