

Probe Covers Lubricated



Clarity Probe Covers are a cutting-edge solution designed for endocavity imaging on ultrasound systems. Crafted from non-latex synthetic polyisoprene, this probe cover ensures a versatile fit across various ultrasound transducers, offering healthcare professionals unmatched flexibility during examinations.

Product Information

- Ultra smooth surface
- Single walled
- Translucent natural colour
- Length : 200 mm \pm 5mm
- Width : 53 mm \pm 2 mm
- Diameter : 34 mm
- Thickness : 0.066 to 0.075mm
- High quality homogenised low allergy rubber latex



Packaging

- 100 probe covers/box
- 20 boxes/carton
- Total pieces/carton - 2000

Carton Size

Box Size

- | | |
|------------------|--------------|
| ■ 37.8x37.8x32cm | ■ 8x6x8cm |
| ■ G.W 7.56kg | ■ Weight 3kg |
| ■ N.W 7.24kg | |

Base Unit



Outer Carton



www.constar.com.au

