

# Image Intensifier Covers



The Task Medical Image Intensifier Cover is designed to be used as a sterile cover for radiographic equipment.

## Product Information

- 90 x 75cm
- Double packed
- Packed into individual pouches
- Latex free
- Single use only
- EO sterilised
- Peel pouch



## Packaging

- 20 pouches/box
- 6 boxes/carton
- Total pieces/carton - 120

## Carton Size

- 55x34x54cm
- G.W 7kg
- N.W 5kg

## Box Size

- 35x20x15cm
- Weight 0.95kg

Base Unit



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

