

Task Patient Gown Half Sleeves - XL



Introducing the SMS Patient Gown with Half Sleeves – designed for optimal patient comfort and medical professional convenience. Crafted from trilaminare non-woven fabric, this gown ensures a perfect balance of opacity and high-level protection.

Product Information

- 10 gowns / bag
- Short sleeves with ties at hip & neck
- SMS - 40GSM
- Hospital Grade
- Single use only
- X-Large - 130x180cm

Packaging

- 10 gowns/bag
- 5 bags/carton
- 50 gowns/carton

Carton Size

- 35x32x33cm
- G.W 5.5kg
- N.W 6kg



CODE
TKGOWSMSSXL

Base Unit



Inner Box



Outer Carton



www.constar.com.au

