

Patient Positioning Pads Solutions



*Enhancing Comfort &
Support in Medical Care*

TASK
MEDICAL

Task Medical Patient Positioning Pad

Designed to optimise patient positioning during female pelvic ultrasounds, particularly on non-tilting couches. Made from high-density foam and encased in hospital-grade vinyl, it offers enhanced comfort and stability with easy-to-clean, durable surfaces for repeated use.

- Breathable fabric
- Antibacterial
- Waterproof and water-repellent zip
- Easy to clean
- Anti-slip fabric




Order Code
TKPOSTV

Streamlined design aids in quicker patient setup and positioning!

Task Medical Foam Positioning Pad 32° Short

This compact 32° foam wedge is crafted to ensure stability during imaging, helping maintain patient positioning for clear and consistent results.




Pack of 2 for an extra on hand

Order Code
TKPOS32

Task Medical Foam Positioning Pad 32° Long

With an extended 32° angle, this positioning pad offers added support along longer sections of the body, ideal for extended imaging procedures.




Pack of 2 for an extra on hand

Order Code
TKPOS32-60

Task Medical Foam Positioning Pad 45° Short

A 45° angled support designed to fit snugly and maintain precise positioning, ensuring reliable imaging outcomes.




Pack of 2 for an extra on hand

Order Code
TKPOS45

Task Medical Foam Positioning Pad 45° Long

Providing a stable 45° incline, this longer pad supports more extensive positioning needs, minimising patient movement for accurate diagnostics.



Pack of 2 for an extra on hand

Order Code
TKPOS45-60



Unit 4, 68 Parramatta Road, Underwood QLD 4119

☎ 1300 778 924 ✉ info@constar.com.au 🌐 www.constar.com.au


Task Medical Foam Positioning Pad Set 8 Pieces

This versatile set of positioning pads includes eight unique shapes designed to support a variety of patient positioning needs during medical imaging or therapeutic procedures. Each pad is made from high-density, radiolucent foam, ensuring durability, patient comfort, and clear imaging.


Features:

- Comprehensive set to accommodate multiple positioning needs
- High-density foam for stability and comfort
- Radiolucent design allows for clear imaging without interference
- Lightweight and easy to arrange for various medical procedures
- Durable and easy to clean, ensuring hygienic and repeated use


**Triangle
45° Pad**



**Finger Step
Pad**



Skull Rest Pad



Circle Pad




**2 x Small
Triangles**



Includes:

- 1 x Triangle 45° Pad - **255 x 180 x 205mm**
- 2 x 15° Wedges - **210 x 47 x 255mm**
- 1 x Circle Pad - **200 x 50mm**
- 1 x Finger Step Pad - **150 x 125 x 57mm**
- 1 x Skull Rest Pad - **255 x 205 x 50mm**
- 2 x Small Triangles - **205 x 75 x 75mm**

2 x 15° Wedges




**Order Code
TKPOSPAD**

Task Medical Uniblock Positioning Device

The Uniblock is a versatile support solution for the head, arms, legs, and feet, designed for physiotherapists, paramedics, and medical professionals. Its lightweight and portable construction allows for easy transport and cleaning, making it an ideal choice for clinical settings. With five different elevations available, it enhances patient positioning while maintaining comfort.

- Hygienic polyurethane with a disinfectable skin
- CE marked for safety compliance
- Supports multiple body parts with adjustable heights
- Lightweight design (300g) for easy handling
- Easily cleaned with medical disinfectant or autoclaved



**Order Code
TKBBLOCK**

Ergonomic design simplifies patient support and optimises comfort during treatment!



Unit 4, 68 Parramatta Road, Underwood QLD 4119

☎ 1300 778 924 ✉ info@constar.com.au 🌐 www.constar.com.au



Get in Touch!




1300 778 924



info@constar.com.au



www.constar.com.au



Scan me!

