

Deluxe Gynaecology Chair - Navy Blue

TASK
MEDICAL

The Task Medical Deluxe Gynaecology Chair in Navy Blue is designed to provide superior comfort and support for both patients and practitioners. Featuring adjustable stirrups, footboards, and an ergonomic reclining mechanism, it ensures optimal positioning for gynaecological examinations and procedures. Its flexible adjustments accommodate various clinical needs, enhancing accessibility and efficiency. With a sturdy build and easy-to-clean upholstery, this chair is a durable and practical choice for medical facilities.

Product Features

- Double gas strut back adjustment (0 to 80°)
- Electric height & seat adjustment (590–1060mm)
- Four retractable medical-grade castors
- Adjustable armrests, stirrups & footboards
- Stainless steel debris tray for hygiene
- Medical-grade PVC vinyl upholstery
- Supports up to 250kg

CODE
TKGCNB



Product Warranty

- Frame: 5 years
- Electrical components: 3 years
- Upholstery: 1 year

Outer Carton



www.constar.com.au

