

EZYM Multipurpose Mobility Bed

TASK
MEDICAL

The Task Medical EzyM Multipurpose Mobility Bed is a manual wind hi-lo patient mobility bed, designed for versatile use in medical practices, aged care facilities, and hospitals. Its lightweight yet durable construction offers easy manoeuvrability and patient comfort.

Product Features

- Manual wind hi-lo adjustment (56 cm – 86 cm)
- Gas spring backrest (tilt 0 – 60°)
- Quick-access drop sides
- Centre-mounted fifth castor for enhanced manoeuvrability
- Four-hook IV pole (adjustable height, detachable)
- Hospital-grade vinyl mattress (60 mm thick)
- Safe working load: 180 kg
- 150 mm quiet-roll castors with front foot brakes
- Two oxygen cylinder holders
- Detachable side rails
- Dimensions: L 197 cm x W 64 cm

CODE
TKEZYM



Outer Carton



Ideal For Use As

- Examination couch
- Procedure table
- Emergency stretcher
- Patient transporter
- Recovery room bed



www.constar.com.au

