

Easy Roller Swimming Pool Wheelchair



The EasyRoller Swimming Pool Wheelchair is completely metal-free, making it ideal for pools, spas, and rehabilitation centres. Its corrosion-resistant design ensures durability and hygiene in wet environments. Developed with NTNU and Eker Design, it meets industry standards, including MRI-safe certification. With an innovative drainage system and ergonomic build, it offers a safe and reliable mobility solution for aquatic use.

Product Features

- 100% metal-free & rust-proof
- Durable materials: rotomoulded PE, nylon, PU tyres
- Water drainage system for quick drying
- Stability-enhancing stainless steel weights
- User-friendly, no loose parts, easy to manoeuvre
- Minimal maintenance, built for institutional use
- Tip-stop safety feature prevents tipping beyond 5°

CODE
WHEEASYRLSW



Product Specifications

- Weight: 17 kg
- Dimensions: 67 x 82.5 x 105 cm (W x L x H)
- Warranty: Body & Footrest – 5 years, Other parts – 1 year



www.constar.com.au

