

## SS Trolley Single Three Drawer with Lock

**TASK**  
MEDICAL

This version of the Task Medical Single Three Drawer Trolley includes a secure lock on the top drawer for added safety. Ideal for environments requiring both storage and security, it combines durable construction with easy mobility.

### Product Features

- Key lock on top drawer for security of contents (spare key included)
- High-grade 304 stainless steel throughout
- Soft-closing magnetic drawers
- German-made ball-bearing drawer runners
- 125mm deep drawers
- 95mm smooth-rolling casters for easy mobility
- Locks on front casters
- Top and bottom side rails on all three sides
- Double safety folds
- Smooth welding finish

### Product Dimensions

- (W) 50cm (D) 50cm (H) 90cm

CODE  
TK3D5OLK



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

