

Clinical Trolley

TASK MEDICAL

The TASK Medical Clinical Trolley is a lightweight, manoeuvrable utility trolley with four large castors and front toe-lock brakes for smooth movement and stability. Built with a heavy-duty aluminium frame and ABS plastic, it offers a durable, ultra-stable workstation suitable for various medical and utility applications.

Product Features

- Lightweight yet sturdy aluminium frame with ABS body
- Top workbench with recessed lip edge and non-slip plastic mat
- Full-width drawer with adjustable dividers and smooth steel runners
- Two under shelves for additional storage
- Three-sided stainless steel handrail for easy manoeuvring
- Two dust bins (grey and yellow flip lids) for waste disposal
- Extra-large castors with two front brakes for stability

CODE
TKTCID



Product Specifications

- Material: Hospital-grade ABS, aluminium frame, stainless steel guard rail
- Drawer: One, with dividers
- Shelves: Two under shelves
- Casters: Roller castors, two with brakes
- Dimensions: 90H x 60W x 50D cm

Outer Carton



www.constar.com.au

