

## Double Step Stool



The Task Medical Double Step Stool is designed for healthcare and home use, helping patients step onto and down from medical examination couches or beds. It's also ideal for reaching shelves or cupboards in slightly higher-than-normal applications, providing a stable alternative to a ladder. With a wider, deeper standing platform and supportive handrail, this stool prioritises stability, making it especially helpful for those with physical challenges or recovering from surgery.

### Product Features

- Double step design for stable access
- One side handrail (fits either side) for added support
- Powder-coated steel frame for strength
- Non-slip wooden steps for secure footing
- Feet caps to protect floors
- Ships in two sections – base and handrail
- Easy to assemble
- Meets OH&S requirements

CODE  
TKDSTEPHR



### Product Specifications

- Step width: 60 cm
- 1st Step depth: 30 cm
- 2nd Step depth: 40 cm
- 1st Step height: 22 cm
- 2nd Step height: 42 cm
- Handrail height: 97 cm (lower), 42 cm (higher)
- Frame diameter: 22 mm

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

