

Emergency Trolley with Locking Bar

TASK
MEDICAL

The Task Medical Emergency Trolley is a robust and efficient solution designed for Australian healthcare facilities. Featuring a 5-drawer configuration, locking bar, defibrillator shelf, and IV pole, it ensures quick and responsive manoeuvrability in emergencies. The latest model includes off-white surface tops, sides, and back, with improved drawer handles.

Product Features

- Secure locking bar and central key locking system
- Durable alloy frame with a washable ABS body
- Five soft-closing drawers with removable dividers and label cards
- Waste bins, needle disposal flask, utility tubs, and retractable shelf
- CPR board, oxygen tank holder, and adjustable defibrillator shelf
- Extra-large noiseless casters with front brakes



Product Specifications

- Colour: Red & White
- Dimensions: 628 x 475 x 950mm
- IV Pole: 304 Stainless Steel, foldable 4 hooks, adjustable height
- Power Board: 240V / 10A capacity, 5m cable

Outer Carton



www.constar.com.au

