

Round Stool with Ring - Navy Blue

TASK
MEDICAL

Featuring an under-seat ring mechanism for height adjustment, this round stool combines simplicity with functionality. The padded PU seat ensures comfort, while the sturdy chrome base and PU castors allow for smooth rolling and swivel motion. Ideal for procedural rooms, labs, or any setting where compact, height-adjustable seating is required.

Product Features

- Mech: Ring mechanism
- Chrome gas lift: 16cm
- Foam thickness: 13mm
- Chrome base 500mm diameter
- Double wheel castors
- Weight capacity: 160kg

Carton Size

- 49x20x46.5cm
- Chair Weight: 6.5kgs
- N.W: 6.00kg
- G.W: 6.70kg

CODE
TKIO24NB



Outer Carton



www.constar.com.au

