

SS Mobile IV Drip Stand 4 Hooks

TASK
MEDICAL

The Task Medical IV Drip Stand Economy is designed for ease of use in healthcare settings. Built with high-quality stainless steel, it ensures durability and stability for medical professionals during patient care. The height-adjustable pole, sturdy base, and smooth mobility make it a reliable and practical solution for IV drip administration.

Product Features

- 4 Stainless Steel Hooks: For secure attachment of IV bags.
- Stainless Steel Pole: Durable and easy to clean.
- 5-Leg Stainless Steel Base: Provides a stable foundation to prevent tipping.
- Height Adjustable: Easy height adjustments with a convenient knob.
- 5 Castors: Ensures smooth and easy mobility around healthcare environments.



Outer Carton



www.constar.com.au

