

## Saddle Stool - Black

**TASK**  
MEDICAL

This ergonomic saddle stool is designed to encourage a healthy seated posture by supporting the natural curvature of the spine. Its high-density foam saddle seat, covered in durable PU, reduces pressure on the lower back and is ideal for professionals working long hours. Fitted with a chrome base and gas lift, this stool provides mobility and comfort in clinical and therapeutic settings.

### Product Features

- Mech: Simple mechanism
- Chrome gas lift: 16cm
- Foam thickness: 13mm
- Chrome base 500mm diameter
- Double wheel castors
- Weight capacity: 160kg

### Carton Size

- 49x20x46.5cm
- Chair Weight: 6.5kgs
- N.W: 5.20kg
- G.W: 5.90kg

CODE  
TK3007-4BLK



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

