

Saddle Stool with Back - Black

TASK
MEDICAL

Combining a supportive saddle seat with an elegant backrest, this stool is built for comfort and control during long procedures. Its high-density foam and PU finish provide durability and ease of cleaning. Ideal for clinical, imaging, and consultation use where posture and movement are key.

Product Features

- Mech: Multi-Fun mechanism
- Chrome gas lift: 16cm
- Foam thickness: 13mm
- Chrome base 600mm diameter
- Double wheel castors
- Weight capacity: 160kg

Carton Size

- 58x25x55cm
- Chair Weight: 8.5kgs
- N.W: 8.90kg
- G.W: 9.70kg

CODE
TKIO37-3BLK



Outer Carton



www.constar.com.au

