

Single Step Stool



The Task Medical Single Step Stool is designed for both healthcare environments and home use, providing assistance for patients when stepping onto and down from medical examination couches or beds. With a wider and deeper platform, this stool offers enhanced stability. The addition of a handrail aids in safe stepping, making it suitable for patients recovering from surgery or those with physical challenges.

Product Features

- Single step design for easy access
- One side handrail for support and safety
- Sturdy tubular powder-coated steel frame for durability
- Non-slip surface wooden step tread for secure footing
- Feet caps to protect floors
- Ships in two sections – base and handrail
- Easy to assemble
- Helps meet workplace OH&S requirements

CODE
TKSSTEPHR



Product Specifications

- Dimensions: 60W x 30D x 29H cm (Step)
- Handrail height: 112 cm
- Tubular frame diameter: 22 mm
- Step width: 60 cm
- Step depth: 30 cm
- Step height: 29 cm

Outer Carton



www.constar.com.au

