

## Step Stool



The Task Medical Step Stool offers a sturdy and reliable solution for reaching higher places. Its chromed steel construction and reinforced rubber tips ensure both durability and safety. Designed for ease of use, it provides stability with a non-skid rubber platform, making it perfect for medical and general environments.

### Product Features

- Chromed steel construction for durability
- 1" steel tubing for strength
- 150 kg capacity for robust support
- Non-skid rubber platform for secure footing
- Reinforced rubber tips for added stability
- Silver vein colour for a sleek appearance
- Heavy-duty for long-lasting use

CODE  
TKSTOIF



### Product Specifications

- Dimensions: 36 x 28 x 22 cm

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

