



“PARR” Automatic Blood Pressure Monitor

Model: Z5



♥ PARR

♥ AFib Atrial Fibrillation Detection (AFib)

The upper chambers of the heart (the atria) beat irregular (quiver) and mostly fast instead of beating effectively to move blood into the ventricles. This condition is associated with a higher risk for heart-related complications. About 10 % - 20 % of patients who suffer from an ischemic stroke also suffer from atrial fibrillation.

♥ PC Premature Contraction Detection (PC)

Extra, abnormal heartbeats generated in abnormal locations of your heart, either in the atria (PAC) or in the ventricles (PVC) or the cardiac conduction nodes (PNC). These extra beats disrupt your regular heart rhythm, sometimes causing palpitations (e.g. skipped beats) in your chest. It may occur singularly or repetitively with various incidence.

♥ PARR Technology

You can prevent the risk of cardiac disorders by early detection. Cardiac disorders need medical attention and thus early diagnosis is of paramount importance. The PARR technology detects arrhythmia during regular blood pressure checks without any additional user skills, user interaction and measurement prolongation.

- Bluetooth® low energy Technology
- Arrhythmia Detection (ARR)
- Atrial Fibrillation Detection (AFib)
- Premature Contraction Detection (PC)
- Cuff Detection
- Movement Detection
- Hypertension Risk Indication
- Average of last 3 readings
- 7 day AM/PM Averaging
- 120 Memories for 2-Users with Date & Time Indication
- Guest Mode
- Large LCD Size : 69*64.4 mm
- USB Type C port (USB Power Connector)
- Rechargeable Lithium-ion Battery
- Lifetime Calibration
- Zipper Bag
- Latex-Free Patented Universal Cone Cuff



Bluetooth® low energy Technology



ARR/AFib



ARR/PC



Cuff Detection



Movement Detection



Hypertension Risk Indication



Average of last 3 readings



7 Day AM/PM Averaging



Guest Mode



120 memories for 2-users with date and time indication



Lifetime Calibration



Type C Port



Universal Cone cuff fits arm circumference: 24~40cm (9.4"~15.7") Other size cuffs available upon request



Important

The device is not intended to replace traditional methods of diagnosis. If arrhythmia (ARR), premature contraction (PC) or atrial fibrillation (AFib) is detected consult a doctor for potential diagnosis or treatment. If user has symptoms of heart attack, stroke or other cardiovascular conditions do not rely on the notification of the device and consult a doctor.

Distributed by:

Silverline Medical Pty Ltd

8A, 380 Pennant Hills Road, Pennant Hills, NSW 2120

1800 959 777

www.rossmax.com.au