

Fingertip Pulse Oximeter

Model: SB220



- Instant readings of SpO₂, Pulse rate
- Shield Design against ambient light
- Biocompatibility protection
- Pulse Strength Indicator
- Visible & audible alert
- Two way, two color OLED display
- String-attached



%SpO₂ & Pulse Rate



Shield Design against Ambient Light



Biocompatibility protection



Pulse Strength Indication



Visible/Audible Alert



OLED Two-way Two-Color Display



String-Attached

Specification

- Measurement range: SpO₂: 35%~99%; PR: 30~250 bpm
- Precision: SpO₂: 70%~99%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3
- Alarm: Default Value for SpO₂
- Power Supply: 2"AAA" Alkaline

Model	Qty per carton	Carton volume
SB220	40 pcs	0.028 cbm/ctn

● Storage temperature: -25°C - 70°C (-13°F - 158°F)
Relative humidity: 15% - 90% (no condensing)
Atmospheric Pressure: 700 to 1060 hPa

Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring