





disposeal+ ACTIVE SAFETY IV CANNULA (CATHETER) with blood control without wings and without injection port



TECHNICAL SPECIFICATION :

| Size | Catheter OD (mm) | Catheter length (mm) | Flow Rate (Water) (ml/min) | Box Qty. (Inner/Outer) |
|------|------------------|----------------------|----------------------------|------------------------|
| 14G | 2.1 | 32 | 270 | 50/500 |
| 16G | 1.8 | 32 | 200 | 50/500 |
| 18G | 1.3 | 32 | 90 | 50/500 |
| 20G | 1.1 | 32 | 60 | 50/500 |
| 22G | 0.9 | 25 | 36 | 50/500 |
| 24G | 0.7 | 19 | 23 | 50/500 |
| 26G | 0.6 | 19 | 15 | 50/500 |

CERTIFICATION :

•CE •ISO 13485:2016