



### FEATURES :

**NEEDLE :** Siliconized, Ultra Sharp, Multi facet, stainless steel back cut, needle help clinician to perform easy and painless venipuncture.

**CATHETER :** Available in FEP, PUR and PTFE material. Siliconized with proprietary technology ensuring smooth surface and minimal discomfort for both patient and clinician.

**FLASHBACK CHAMBER :** Large FB chamber allow quick visualisation for successful venipuncture.

**NEEDLE COVER :** Transparent needle cover to prevent accidental damage to needle and catheter.

**FLOW CONTROL HUB :** The flow control hub has hydrophobic filter for allowing air to pass and don't allow blood to ooze out.

**CAP :** Colour Coded cap for gauge identification.

**OPTION AVAILABLE :** FEP/PTFE/PUR Catheter | Hydrophobic filter | Quick flash back needle  
Packing colour coded M.K.P. | Sterile/non-pyrogenic/Single use.

### TECHNICAL SPECIFICATION :

**CERTIFICATION :** •CE •ISO 13485:2016

Size	Catheter OD (mm)	Catheter length (mm)	Flow Rate (Water) (ml/min)	Box Qty. (Inner/Outer)
14G	2.1	45	270	50/500
16G	1.8	45	200	50/500
17G	1.5	45	125	50/500
18G	1.3	45	90	50/500
20G	1.1	32	60	50/500
22G	0.9	25	36	50/500
24G	0.7	19	23	50/500
26G	0.6	19	15	50/500