



CLARITY

LEADING THE WAY IN ULTRASOUND TECHNOLOGY

Our extensive portfolio of ultrasound products is used and trusted by Sonographers throughout Australia. Experience superior quality and performance.

CONS  AR



1300 778 924

Find out more at www.constar.com.au

Clarity Ultrasound Transmission Gel

**Advanced
Formula!**

Our Clarity Ultrasound Transmission Gel is designed for superior acoustic performance. **Available in both 5L and 250ml bottles, in blue or clear,** it ensures clear, uninterrupted ultrasonography and Doppler applications.

- Suitable for all ultrasound frequencies
- Hypoallergenic, CIT/MIT, and paraben-free
- High viscosity, vacuum-treated
- Water-soluble, oil-free, odourless



Blue

Clear

**250ML
12 Pack**

CLUSG250BL



CLUSG250CL

**Fresh Look
Same Quality**



5 Litre

**Blue
CLUSG5BL**

**Clear
CLUSG5CL**

Clarity Gel Warmers

**Maintain Constant
Temperature**

Enhance patient comfort with Clarity Ultrasound Gel Warmers. **Available in single and dual bottle models,** these warmers feature adjustable temperature control and a sleek design, making them an essential addition to any ultrasound setup.

- LED display, adjustable temperature (34-40°C)
- Programmable timer, energy-efficient
- Compact design

**Perfect Transmission
Of Soundwaves**



**Single
CLGW01**



**Double
CLGW02**



1300 778 924

Find out more at www.constar.com.au

Clarity Probe Covers

**One Size
Fits Most**

Clarity Probe Covers provide reliable protection for your ultrasound probes. **Available in lubricated, non-lubricated, and latex-free options**, these covers ensure hygienic and efficient examinations.

- Single walled
- Ultra smooth surface
- Rolled for convenience



**Lubricated
CLPCVL**



**Non-Lubricated
CLPCV**



**Latex Free
CLPCVLF**



Clarity Freezer Bags

Latex-Free Bands

Our Clarity Sterile Freezer Bags are designed for single-use and come with two latex-free elastic bands. These bags offer a sterile and convenient solution for your medical storage needs.

- 30 x 45cm
- Includes 2 latex-free elastic bands
- Individually packed



30 x 45cm

**Individually
Packed**



**Order Code
CLFBW**



1300 778 924

Email us at info@constar.com.au



Enhance Your Ultrasound Experience with These Essential Products

Discover our carefully curated selection of complementary products designed to enhance the efficiency and effectiveness of your ultrasound procedures.

Positioning Pad

A high-density foam wedge designed for optimal patient positioning during female pelvic ultrasound. It features a durable, hospital-grade vinyl cover that is easy to clean and provides a stable, comfortable surface on non-tilting couches.

- Rigid for support
- Soft for comfort
- Easy to clean

150 x 550 x 500mm

Durable & Hygienic

TASK
MEDICAL

Order Code
TKPOSTV



Positioning Pad 45°

Available in a convenient pack of 2, perfect for always having a backup ready when you need it.


Enhanced Stability

255 x 180
x 205mm

258 x 180
x 60mm



Short 45°
TKPOS45




Long 45°
TKPOS45-60

Positioning Pad Kit

A versatile set of 8 pads, each tailored for specific positioning needs, enhancing imaging accuracy.

Precision Positioning



Order Code
TKPOSPAD



1300 778 924

Find out more at www.constar.com.au