

Cuff for Blood Pressure Monitor



Patented Cone Cuff



Position Guideline



Latex Free



Easy to Apply



Clear Position Guideline



Cone Cuff Patent Nos :
ZL 2003 2 0100815.4
Nr. 20 2004 014 458.1
M244879
Other patents pending



Digital Korotkoff
Cone Cuff



SMALL

18 - 26 cm

7.0"-10.3"



MEDIUM

24 - 36 cm

9.5"-14.2"



LARGE

34 - 46 cm

13.4"-18.1"



UNIVERSAL

24 - 40 cm

9.4"-15.7"

Model	Qty per carton
Cone Cuff / Small	48 pcs
Cone Cuff / Medium	
Cone Cuff / Large	
DK Cuff / Small	48 pcs
DK Cuff / Medium	
DK Cuff / Large	
Universal Cuff	48 pcs
Carton volume	0.088 cbm/ctn

Storage and Transportation Condition
-10°C~60°C(14°F~140°F)
10%~90% RH; 700~1060 hPa

Distributed by:
CONSTAR X PTY LTD
Unit 4, 68 Parramatta Road, Underwood, QLD 4119,
Australia. Ph: 1300 778 924
www.constar.com.au