



“XL” Deluxe Automatic Wrist Blood Pressure Monitor

Model: BQ705



- Jumbo size “XL” digit
- 120 memories for 2-users with date and time indication
- Guest Mode
- Movement Detection
- Smart sense technology
- Hypertension Risk Indication
- Irregular Heartbeat Detection (IHB)
- Average of last 3 readings
- One-touch automatic operation
- Carrying case
- PC Link (Optional)
- Lifetime calibration



Jumbo size
“XL digit”



120 memories
for 2-users



Guest Mode



Movement
Detection



Smart Sense
Technology



Hypertension
Risk Indication



Irregular Heartbeat
Detection (IHB)



Average of
last 3 readings



AAMI

This blood pressure monitor meets the AAMI accuracy standard requirements (AAMI/ANSI/ISO 81060-2)



Cuff fits wrist
circumference: 13.5~22cm
(5.3~8.7")

Smart Sense Technology

Smart Sense Technology utilizes an advanced air control system. The electronic deflation valve precisely controls the deflation rate and allows the monitor to detect and analyze pulse signals, adjusting cuff pressure during deflation. This significantly ensures more precise analysis, while reaching greater accuracy.

Model	Qty per carton	Carton volume
BQ705	36 pcs	0.046 cbm/ctn

Storage and Transportation Condition
-10°C~60°C(14°F~140°F) 10%~90%
RH, 700~1060 hPa

Distributed by:

CONSTAR X PTY LTD

Unit 4, 68 Parramatta Road, Underwood, QLD 4119,
Australia. Ph: 1300 778 924

www.constar.com.au