

Fingertip Pulse Oximeter

Model: SB100



- %SpO2 & Pulse Rate - Big LED screen display for easy reading
- Shielded Design blocks Ambient Light
- Biocompatibility & Anti-Allergic Design
- One- piece face cover design & One-touch ON/OFF key
- Pulse Strength Indicator
- Lanyard-attached - Economical & portable health management tool



%SpO2 & Pulse Rate



Large Font Size



Shielded Design blocks Ambient Light



Biocompatibility & Anti-Allergic Design



Pulse Strength Indication



String-Attached



Specification

Measurement Range | SpO2: 35%~99%; PR: 30~250 bpm

Precision | SpO2: 70%~99%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3

Power Supply | 2 x "AAA"Alkaline

Model	Qty per carton	Carton volume
SB100	40 pcs	0.028 cbm/ctn

Storage temperature: -25°C - 70°C (-13°F - 158°F)

Relative humidity: 15% - 90% (no condensing)

Atmospheric Pressure: 700 to 1060 hPa

Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

Distributed by:

CONSTAR X PTY LTD

Unit 4, 68 Parramatta Road, Underwood, QLD 4119, Australia. Ph: 1300 778 924

www.constar.com.au