

## Fingertip Pulse Oximeter

Model: SB220



- Instant readings of SpO<sub>2</sub>, Pulse rate
- Shield Design against ambient light
- Biocompatibility protection
- Pulse Strength Indicator
- Visible & audible alert
- Two way, two color OLED display
- String-attached



%SpO<sub>2</sub> & Pulse Rate



Shield Design against Ambient Light



Biocompatibility protection



Pulse Strength Indication



Visible/Audible Alert



OLED Two-way Two-Color Display



String-Attached

### Specification

- Measurement range: SpO<sub>2</sub>: 35%~99%; PR: 30~250 bpm
- Precision: SpO<sub>2</sub>: 70%~99%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3
- Alarm: Default Value for SpO<sub>2</sub>
- Power Supply: 2"AAA"Alkaline

Model	Qty per carton	Carton volume
SB220	40 pcs	0.028 cbm/ctn

Storage temperature: -25°C - 70°C (-13°F - 158°F)  
Relative humidity: 15% - 90% (no condensing)  
Atmospheric Pressure: 700 to 1060 hPa

### Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

Distributed by:

**CONSTAR X PTY LTD**

Unit 4, 68 Parramatta Road, Underwood, QLD 4119, Australia. Ph: 1300 778 924

[www.constar.com.au](http://www.constar.com.au)