

rossmax australia

Prevention

Monitoring

Therapy

ROSSMAX DISTRIBUTOR CATALOGUE 2024.3

CONStar



| 90
COUNTRIES



| 85
PRODUCTS



| 51
PATENTS



ACT
ARTERY CHECK
TECHNOLOGY

PARR
PULSE ARRHYTHMIA
TECHNOLOGY

www.rossmax.com.au



HEALTHSTYLE APP
NOW AVAILABLE



Unit 4, 68 Parramatta Road, Underwood QLD 4119
☎ 1300 778 924 ✉ info@constar.com.au 🌐 www.constar.com.au



rossmax
health station



BLOOD PRESSURE
MONITORING



ARTERY CONDITION
CHECK



SPO2
MONITORING



WEIGHT & BODY FAT
MANAGEMENT



TEMPERATURE
MEASUREMENT



5 Professional Health Checks in 1

- ➕ **Enhance** your health professional services with innovative & unrivalled technology
- ➕ **Pioneer** effective screening for major chronic diseases
- ➕ **Prevent** stroke & heart attacks by early detection of risk factors
- ➕ **Boost** customer retention and loyalty with personalised health checks
- ➕ **Endorse** home monitoring with unique ACT & PARR technologies and help save lives
- ➕ **Journey** with Australians to a better Healthstyle®

INCLUDED CONNECTED DEVICES



RMX5BT PARR
Blood Pressure Monitor



RMSB210BT ACT
Finger Pulse Oximeter



RMWF262BT
Digital Bluetooth Scale



RMHC700BT
Telephoto Thermometer



Unit 4, 68 Parramatta Road, Underwood QLD 4119
☎ 1300 778 924 ✉ info@constar.com.au 🌐 www.constar.com.au



rossmax

health station



Bluetooth

Health Station - STM300F Desktop PRO | 15.6"



ARR/AFib



ARR/PC



ARTERY
CONDITION



- Touch screen
- Printer
- Health check stand
- Camera
- NFC reader
- Smart card reader
- QR code / barcode scanner
- Full set of connected devices



Health status receipt includes:

- + Health check measurement results
- + QR Code for scanning results to a phone and integrates to the Rossmax Healthstyle App.
- + Optional discount coupon for Rossmax in-store products

ENQUIRE FOR MORE INFO & PRICING



Bluetooth

Health Station - STM200M Desktop MINI | 10.4"



ARR/AFib



ARR/PC



ARTERY
CONDITION



- Touch screen
- Health check stand
- QR Code for scanning results to a phone and integrates to the Rossmax Healthstyle App
- Option to email results to customer
- Full set of connected devices



Unit 4, 68 Parramatta Road, Underwood QLD 4119

1300 778 924 info@constar.com.au www.constar.com.au



Healthstyle App & Connected Devices

5 in 1

Manage your **Blood Pressure, Artery Condition, Oxygen Saturation, Weight & Temperature** all in one App.

Caring Friends

Allows you to share your measurements with Loved Ones, Caregivers & Healthcare Providers.

Health Data

Allows you to import in-store Health Station results to your phone.

 **BLOOD PRESSURE MONITORING**

 **ARTERY CONDITION CHECK**

 **SPO2 MONITORING**

 **WEIGHT & BODY FAT MANAGEMENT**

 **TEMPERATURE MEASUREMENT**



RMX5BT PARR
Blood Pressure Monitor



RMHC700BT
Telephoto Thermometer



RMWF262BT
Digital Bluetooth Scale



RMSB210BT ACT
Finger Pulse Oximeter



rossmax

monitoring



BP Monitor Upper Arm PARR PRO Professional – AC1000F



early stage
RISK SCREENING
with ONE single
BLOOD PRESSURE
measurement



Trolley
(Optional)



300,000
Measurements
Guaranteed!



ARR/AFib



ARR/PC



ARR/TACH



ARR/BRAD

♥ PARRPRO PROFESSIONAL



Atrial Fibrillation detection (AFib)

The upper chambers of the heart (atria) beat irregularly (quiver) and mostly fast instead of beating effectively to move blood into the ventricles. This condition is associated with a higher risk for cardiac blood clots, stroke, heart failure and other heart-related complications. About 10% - 20% of patients who suffer from an ischemic stroke also suffer from atrial fibrillation.



Premature Contraction detection (PC)

Extra, abnormal heartbeats generated in abnormal locations of your heart, either in the atria (PAC) or in the ventricles (PVC) or the cardiac conduction nodes (PNC). These extra beats disrupt your regular heart rhythm, sometimes causing palpitations (e.g. skipped beats) in your chest. It may occur singularly or repetitively with varying incidence. If not stress-related, they are markers for a multitude of cardiac disorders including elevated ischemic stroke risk with PC.



Tachycardia detection (TACH)

Fast heart rate of more than 100 beats per minute (BPM) in adults. Indicator for cardiac and related disorders (e.g. coronary disease, hyperthyroidism, fever) unless the tachycardia has resulted from physical or mental stress.



Bradycardia detection (BRAD)

Slow heart rate of less than 55 beats per minute in adults. Bradycardia may be related to a multitude of cardiac and related disorders (e.g. electrolyte imbalance, medication, coronary and/or valvular disease) unless caused by endurance training heart adaptation.



Non type – specific Arrhythmia (ARR)

• 300,000 Measurements guaranteed! Ideal for use in hospitals, clinics and other medical facilities

• Arrhythmia detection (ARR)

• Atrial Fibrillation detection (AFib)

• Premature Contraction detection (PC)

• Tachycardia detection (TACH)

• Bradycardia detection (BRAD)

• Lithium ion battery pack lasts up to 600 measurements when fully charged

• Professional valve for varying deflation settings 2.5/4.5/6.5 (mm/Hg)

• Cuff detection

• Movement detection

• Automatic & auscultatory measurements

• 7.5V Adaptor included

• Anti-splash protective cover

• S, M, L cuffs included

• Lifetime calibration



Unit 4, 68 Parramatta Road, Underwood QLD 4119

1300 778 924 info@constar.com.au www.constar.com.au



rossmax
monitoring



BP Monitor Upper Arm PARR PRO Bluetooth – X9BT



Bluetooth



ARR/AFib



ARR/PC



ARR/TACH



Cuff Detection

♥ PARRPRO



Atrial Fibrillation detection (AFib)

The upper chambers of the heart (the atria) beat irregularly (quiver) and mostly fast instead of beating effectively to move blood into the ventricles. This condition is associated with a higher risk for cardiac blood clots, stroke, heart failure and other heart-related complications. About 10% - 20% of patients who suffer from an ischemic stroke also suffer from atrial fibrillation.



Premature Contraction detection (PC)

Extra, abnormal heartbeats generated in abnormal locations of your heart, either in the atria (PAC) or in the ventricles (PVC) or the cardiac conduction nodes (PNC). These extra beats disrupt your regular heart rhythm, sometimes causing palpitations (e.g. skipped beats) in your chest. It may occur singularly or repetitively with varying incidence. If not stress-related, they are markers for a multitude of cardiac disorders including elevated ischemic stroke risk with PC.



Tachycardia detection (TACH)

Fast heart rate of more than 100 beats per minute (BPM) in adults. Indicator for cardiac and related disorders (e.g. coronary disease, hyperthyroidism, fever) unless the tachycardia has resulted from physical or mental stress.



Bradycardia detection (BRAD)

Slow heart rate of less than 55 beats per minute in adults. Bradycardia may be related to a multitude of cardiac and related disorders (e.g. electrolyte imbalance, medication, coronary and/or valvular disease) unless caused by endurance training heart adaptation.



Non type – specific Arrhythmia (ARR)

- Bluetooth® low energy technology
- Arrhythmia detection (ARR)
- Atrial Fibrillation detection (AFib)
- Premature Contraction detection (PC)
- Tachycardia detection (TACH)
- Bradycardia detection (BRAD)
- Cuff detection

- Both automatic and auscultatory measurement available
- Movement detection
- Hypertension risk indication
- 90 Memories with date & time
- S, M, L cuffs included
- Cuff holder design
- Jumbo size XL digit (LCD size: 99 x 115mm)
- Backlight LCD
- Adaptor socket
- Lifetime calibration
- PC link (optional)
- Zipper case



Unit 4, 68 Parramatta Road, Underwood QLD 4119
1300 778 924 info@constar.com.au www.constar.com.au



rossmax

monitoring



BP Monitor Upper Arm PARR Bluetooth - X5BT



Bluetooth



7 Day Averaging



ARR/AFib



ARR/PC

- Healthstyle App
- Bluetooth® low energy technology
- 7 Day am/pm averaging
- Arrhythmia detection (ARR)
- Atrial Fibrillation detection (AFib)
- Premature Contraction detection (PC)
- Cuff detection
- Movement detection
- Hypertension risk indication
- Average of last 3 readings
- 120 Memories for 2-users with date & time indication
- Guest mode
- Large LCD size: 69 x 78mm
- Lifetime calibration
- Zipper case



BP Monitor Upper Arm Deluxe Bluetooth - X3BT



Bluetooth



Fuzzy Technology



Irregular Heartbeat Detection (IHB)



Hypertension Risk Indication

- Healthstyle App
- Bluetooth® low energy technology
- Real Fuzzy Technology
- Irregular Heartbeat detection (IHB)
- Hypertension risk indication
- Average of last 3 readings
- 120 Memories for 2-users with date & time indication
- Cuff detection
- Movement detection
- Latex-free patented universal cone cuff
- One-touch automatic operation
- Cuff holder design
- Large LCD size: 57 x 73mm
- Dual power ready
- Lifetime calibration
- PC link (optional)
- Carry case



Unit 4, 68 Parramatta Road, Underwood QLD 4119

1300 778 924 info@constar.com.au www.constar.com.au



rossmax
monitoring



BP Monitor Upper Arm Deluxe XL - MW701F



Cuff Detection (CD)



Detects and indicates whether the cuff is wrapped around the arm properly

Movement Detection



Movement detection helps remind the user to remain still and indicates any movement during measurement. The specified icon appears once movement has been detected during and after each measurement.



Cuff Detection



Movement Detection



Irregular Heartbeat Detection (IHB)



Hypertension Risk Indication

Real **f**uzzy Technology

Real Fuzzy Technology determines ideal cuff pressure based on one's systolic blood pressure and arm size. Users no longer need to pre-set the inflation level before measurement. The technology eliminates arm discomfort through unnecessary high inflation settings and also prevents inaccurate readings caused by erroneous cuff inflation levels.

- Cuff detection
- Movement detection
- Irregular Heartbeat detection (IHB)
- Hypertension risk indication

- Average of last 3 readings
- 7 day AM/PM averaging
- 120 Memories for 2-users with date & time indication
- Guest mode
- Latex-free patented universal cone cuff
- Large LCD size: 101 x 118mm
- Dual power ready
- Lifetime calibration
- PC link (optional)
- Carry case



Unit 4, 68 Parramatta Road, Underwood QLD 4119
1300 778 924 info@constar.com.au www.constar.com.au



rossmax

monitoring



BP Monitor Upper Arm PARR Lithium-Ion Bluetooth – Z5BT



Bluetooth

AUSTRALIA'S
FIRST BPM WITH
RECHARGEABLE
LITHIUM-ION
BATTERY

Provides
600 readings
from one
charge!



Type C Port



Rechargeable
Lithium-ion battery



ARR/AFib



ARR/PC



7 Day
Averaging

- USB powered (Type C)
- Rechargeable Lithium-ion battery
- Bluetooth® low energy technology
- Arrhythmia detection (ARR)
- Atrial Fibrillation detection (AFib)
- Premature Contraction detection (PC)
- 7 Day AM/PM averaging
- Cuff detection
- Movement detection
- Hypertension risk indication
- Average of last 3 readings
- 120 Memories for 2-users
- Guest mode
- Large LCD size: 69 x 64.4mm
- Lifetime calibration
- Latex-free patented universal cone cuff
- Zipper bag



BP Monitor Upper Arm USB Powered - Z1



Type C Port



Cuff Detection



Movement
Detection



Hypertension
Risk Indication



Irregular Heartbeat
Detection (IHB)

- USB powered (Type C)
- Cuff detection
- Movement detection
- Hypertension Risk Indication (HRI)
- Irregular Heartbeat detection (IHB)
- Real Fuzzy technology
- Average of last 3 readings
- 90 Memories
- One-touch automatic operation
- Large LCD size: 69 x 64.4mm
- Lifetime calibration
- Latex-free patented universal cone cuff
- Carry case



Unit 4, 68 Parramatta Road, Underwood QLD 4119
1300 778 924 info@constar.com.au www.constar.com.au



rossmax monitoring




BP Monitor Upper Arm Deluxe - X3



Fuzzy Technology



Irregular Heartbeat Detection



Average of last 3 readings



120 Memories for 2-users with date and time indication

- Real Fuzzy Technology
- Irregular Heartbeat detection (IHB)
- Average of last 3 readings
- 120 Memories for 2-users with date and time indication
- Cuff detection
- Movement detection
- Hypertension risk indication
- Latex-free patented universal cone cuff
- One-touch automatic operation
- Cuff holder design
- Large LCD size: 57 x 73mm
- Dual power ready
- Lifetime calibration
- Carry case



BP Monitor Wrist Deluxe - BQ705



Smart Sense Technology



120 memories for 2-users



Hypertension Risk Indication



Average of last 3 readings

- Smart Sense Technology
- 120 Memories for 2-users with date and time indication
- Hypertension risk indication
- Average of last 3 readings
- Movement detection
- Jumbo size screen
- Irregular Heartbeat detection (IHB)
- One-touch automatic operation
- PC link (Optional)
- Lifetime calibration
- Carry case





rossmax
monitoring

Blood Pressure Monitor - Cuffs

**Compatible
with most
brands**



Available In: S, M, L & Universal

The ergonomically designed Rossmax Cone Cuff and D-Ring makes it easier for an individual to apply by themselves without the cuff slipping or sliding.



SMALL	MEDIUM	LARGE	UNIVERSAL
18 - 26 cm	24 - 36 cm	34 - 46 cm	24 - 40 cm
7.0" - 10.3"	9.5" - 14.2"	13.4" - 18.1"	9.4" - 15.7"



CODE:

PRODUCT:

RMX9BT	BP Monitor Upper Arm - PARR PRO Bluetooth
RMX5BT	BP Monitor Upper Arm - PARR Bluetooth
RMX3BT	BP Monitor Upper Arm - Delux Bluetooth
RMX3	BP Monitor Upper Arm - Delux
RMZ5BT	BP Monitor Upper Arm - PARR Lithium-Ion Bluetooth
RMZ1	BP Monitor Upper Arm - USB Powered
RMBQ705	BP Monitor Wrist - Deluxe
RMMW701F	BP Monitor Upper Arm - Deluxe XL
RMBPCF001	BP Cone Cuff Large (34-46cm)
RMBPCF002	BP Cone Cuff Medium (24-36cm)
RMBPCF003	BP Cone Cuff Small (18-26cm)
RMBPCF004	BP Cone Cuff Universal (24-40cm)



Unit 4, 68 Parramatta Road, Underwood QLD 4119
 1300 778 924 info@constar.com.au www.constar.com.au



Aneroid Sphygmomanometer - GB102



- White enamel 300 mmHg no-pin stop manometer
- Chrome-plate brass air release valve
- Latex free Cuff
- Single head Stethoscope
- Zipper bag



Palm Type Sphygmomanometer - GD102



- Superior performance and quality
- Large black dial face, 65mm easy-to-read
- Single-hand operation
- 300mmHg no-pin stop manometer
- Latex free cuff
- Single head stethoscope
- Zipper bag



CODE:

RMGB102
RMGD102

PRODUCT:


Sphygmomanometer Aneroid
Sphygmomanometer Palm Type



rossmax
monitoring



Infrared Non-Contact Telephoto Thermometer - HC700BT



Bluetooth

1"




- Healthstyle App
- Bluetooth® low energy technology
- Telephoto temperature measurement distance within 10cm (3.94 inch)
- 1 Second measurement
- Temple temperature
- Object temperature
- Auditory fever alarm
- 9 Memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction

Infrared Non-Contact Temple Thermometer - HA500

Quiet



1"



- Measurement distance within 5cm
- 1 Second measurement
- Temple temperature
- Object temperature
- Fever alarm
- 9 Memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator



Unit 4, 68 Parramatta Road, Underwood QLD 4119

1300 778 924 info@constar.com.au www.constar.com.au



rossmax
monitoring



Infrared Ear & Object Thermometer - DET105

with Probe Covers



Storage Case



1^{sec}



- 1 Second measurement
- Ear & object temperature
- 15 Probe covers included
- Storage case included
- Fever alarm
- 10 Reading memory
- °C / °F switchable
- Backlight
- Low battery indicator
- Auto shut-off



CODE:

PRODUCT:

RMHC700BT

Infrared Non-Contact Telephoto Thermometer

RMHA500

Infrared Non-Contact Temple Thermometer

RMHA500Q

Infrared Non-Contact Temple Thermometer

RMDET105

Infrared Ear & Object Thermometer

RMDET105PC


Thermometer Probe Covers 10's

Flexible Tip Thermometer - T4333



10^{sec}

- 10 Second measurement
- Flexible tip
- Fever alarm
- Waterproof
- Auto shut-off
- °C / °F switchable
- Low battery indicator



CODE:

PRODUCT:

RMT4333

Flexible Tip Thermometer 10sec



Unit 4, 68 Parramatta Road, Underwood QLD 4119


1300 778 924




info@constar.com.au



www.constar.com.au



ACT ARTERY CHECK TECHNOLOGY



PULSE OXIMETER with ACT detects

- **ARTERY CONDITION**
- **SpO₂**
- **PULSE RATE**

AT YOUR FINGERTIP

- **STATE OF THE ART**
Uses Artery Check Technology (ACT)
- **EFFORTLESS**
Detection at your fingertip
- **FAST**
Instant readings of SpO₂, pulse rate and artery condition
- **EASY TO UNDERSTAND**
Comprehensive 6-level artery condition indicator

Fingertip Pulse Oximeter - SB210BT

Artery Condition Detection

Artery condition is associated with the potential presence of the following diseases: Arteriosclerosis, peripheral circulation disorder and a wide range of cardiovascular diseases. Monitoring your arterial constriction condition with the Rossmax ACT-embedded Technology allows for an early risk assessment of wide-spread clinical cardiovascular disorders. With this non-invasive technique, both the progression and the actual status of pathological arterial perfusion is recognised. In addition, certain Rossmax blood pressure monitors offer PARR Technology, which is a world first stroke screening method. PARR and ACT technologies combine offering the total monitoring solution to your health by evaluating the most severe health risk factors.

WORLD'S FIRST
ARTERY CONDITION CHECK
WITH SpO₂ & PULSE RATE
PULSE OXIMETER



Bluetooth



ARTERY
CONDITION



SpO₂
PR



- **Healthstyle App**
- **Bluetooth® low energy technology**
- **Artery Check Technology (ACT) embedded**
- **One minute reading for SpO₂, pulse rate & artery condition**
- **Pulse strength indicator**
- Shielded design blocks ambient light
- Two colour OLED display
- BLE 4.0 data transfer real-time



rossmax
monitoring



- Artery Check Technology (ACT) embedded
- Pulse strength indicator
- Instant readings of SpO₂, pulse rate and artery condition in 1 minute
- Shielded design blocks ambient light
- Biocompatibility & anti-allergic design
- Visible & audible alarm
- Two way, two color OLED display
- Cord attached

Fingertip Pulse Oximeter with ACT - SB200



CODE:	PRODUCT:
SB210BT	Fingertip Pulse Oximeter ACT Bluetooth
SB200	Fingertip Pulse Oximeter
SB100	Fingertip Pulse Oximeter

Fingertip Pulse Oximeter - SB100



- SpO₂ and Pulse Rate (PR) readings
- Big LED screen display for easy reading
- Shielded design blocks ambient light
- Pulse strength indicator
- Biocompatibility & hypo-allergenic design
- One-touch ON/OFF key
- Cord attachment





rossmax
monitoring



Handheld Pulse Oximeter - SA300BT



- Bluetooth® low energy technology
- Artery Check Technology (ACT) embedded
- One minute reading for SpO₂, pulse rate & artery condition
- Pulse strength indicator
- Large 2.4" colour TFT display
- Personal SpO₂, heart rate and threshold settings
- Shielded design blocks ambient light
- Biocompatibility protection
- Alarm On/Off
- Silent mode
- Up to 999 memories

Pulse Oximeter Probes



Adult



Neo-Natal



Paediatric



Tube

CODE:	PRODUCT:
RMSA300BT	Handheld Pulse Oximeter - Bluetooth
RMSA120	Handheld Pulse Oximeter
RMPA100	Pulse Oximeter Probe - Adult
RMPB100	Pulse Oximeter Probe - Neo-Natal
RMPC100	Pulse Oximeter Probe - Paediatric
RMPF100	Pulse Oximeter Probe - Tube

- SpO₂ & pulse rate
- Personal SpO₂, heart rate and threshold settings
- Pulse strength indicator
- Shielded design blocks ambient light
- Backlight LCD display for easy reading
- Biocompatibility protection
- Alarm On/Off
- Silent mode
- Up to 288 memories



Handheld Pulse Oximeter - SA120



Unit 4, 68 Parramatta Road, Underwood QLD 4119
1300 778 924 info@constar.com.au www.constar.com.au



rossmax

therapy



Adult with
Cap Only



Infant with
Small Mask



Child with
Medium Mask



Adult with
Large Mask

Aero Spacer - AS175

Valved holding chamber



CODE:

RMA175-C
RMA175-S
RMA175-M
RMA175-L
RMPF120C

PRODUCT:

Aero Spacer - Cap Only
Aero Spacer - S Size Mask
Aero Spacer - M Size Mask
Aero Spacer - L Size Mask
Peak Flow Meter - Child

- **High efficiency one way valve for easy inhalation and exhalation**
- **Universal backpiece, can be used with all commonly prescribed Metered Dose Inhalers (MDI)**
- **Anti-static spacer improves the delivery of MDI by inhibiting the adherence of drug to spacer wall**
- **Soft and flexible silicone mask provides comfortable touch and reliable seal**
- This product contains no latex and is BPA free
- Infant, Child and Adult face masks to cater to the different needs
- Infant (S) mask: 0 - 1.5 yrs
- Child (M) mask: 1 - 5 yrs
- Adult (L) mask or mouthpiece: 5 yrs +



Unit 4, 68 Parramatta Road, Underwood QLD 4119


☎ 1300 778 924 ✉ info@constar.com.au 🌐 www.constar.com.au



rossmax
therapy



Piston Nebuliser - NA100



Efficient Treatment



Patented Bottle



Built-in Compartment for Accessories



Adult and Child Masks

- MMAD \leq 2.2 μ m; Fine Particle Dose (FPD): 80-85%
- Consistent fine particle size for efficient respiratory treatment
- Patented valve adjustable nebuliser bottle (VAT)
- Includes mouthpiece and masks for adult & child

NEBULISER KIT INCLUDES:

- VAT bottle (model N1)
- Mouthpiece
- Filter (5 pieces)
- Air tube (150cm)
- Masks for adult & child

Piston Nebuliser - NI60Q

Qutie



CODE:	PRODUCT:
RMNA100	Piston Nebuliser
RMNI60Q	Piston Nebuliser - Qutie



NEBULISER KIT INCLUDES:

- VAT bottle (model N1)
- Mouthpiece
- Filter (5 pieces)
- Air tube (150cm)
- Mask for child



Unit 4, 68 Parramatta Road, Underwood QLD 4119
☎ 1300 778 924 ✉ info@constar.com.au 🌐 www.constar.com.au




rossmax
therapy

Larger particle size to treat the respiratory tract

Rhinitis
Sinusitis
Pharyngitis
Tonsillitis
Laryngitis
Tracheitis
Tracheobronchitis

COMPLETE
Respiratory Solution

Asthma
Bronchitis
Bronchiolitis
Bronchiectasis
Bronchopneumonia



5 - 10 μm



$\leq 2.2 \mu\text{m}$

Small particle size to reach the deepest lung periphery



Nebuliser - Accessories

VAT Type Accessory Pack

Compatible with most nebuliser models

NEBULISER KIT INCLUDES:

- VAT bottle (model N1)
- Mouthpiece
- Filter (5 pieces)
- Air tube (150cm)
- Masks for adult & child

CODE:

RMNBAC03
RMNBAC04
RMNBAC05
RMNBAC06
RMNBAC02

PRODUCT:

Nebuliser Adult Mask Only
Nebuliser Child Mask Only
Nebuliser Tube Only
Nebuliser Air Filters Only - 5pcs
Nebuliser VAT Type Mask Accessory Pack



Unit 4, 68 Parramatta Road, Underwood QLD 4119

1300 778 924

info@constar.com.au



www.constar.com.au



rossmax
therapy



Portable Suction Unit - V7

Ideal for tracheostomised patients, minor surgical applications and post operative therapy in hospital or home care.

Triple Protection

- 1 Overflow valve 2 Antibacterial filter 3 Tamper-proof design



Order Code: RMV7



Order Code: RMV3

Dual Power

Rechargeable Lithium-ion battery

Dual Power

The V7 consists of a high capacity, rechargeable lithium-ion battery. The high performance pump delivers 600mmHg vacuum and 30L/min flow rate and 25L/min under power-saving mode.

Triple Protection

Overflow valve plus anti-bacterial filter and tamper-proof design prevents back flow of fluid and airborne contaminants from entering the pump.

Convenient

High performance light-weight lithium-ion battery offers approximately 45 minutes of run time under normal mode and up to approx. 60 minutes under power saving mode. The V7 suction unit is mobile and the LED battery indicator shows the remaining power level. There is a storage compartment for all tubes and accessories.

Ergonomic

Comprehensive, ergonomically designed dashboard with button to open the lid.

Portable Suction Unit - V3

Comfortable

High performance pump delivers 550mmHg vacuum and 16L/min flow rate.

Triple Protection

Overflow valve plus anti-bacterial filter and tamper-proof design prevents back flow of fluid and airborne contaminants from entering the pump.

Convenient

There is a storage compartment for all tubes and accessories.

Ergonomic

Comprehensive, ergonomically designed dashboard with button to open the lid.



Unit 4, 68 Parramatta Road, Underwood QLD 4119

1300 778 924 info@constar.com.au www.constar.com.au




Anti Decubitus Mattress Cell Type - AM40



Adjustable hangers to suit any kind of hospital bed frame

Order Code: RMAM40

Accessories




- **2 in 1 alternating therapy**
- **Comfortable nylon cell mattress**
- **Weight indicator with pressure adjustment settings**
- **Low operating sound**
- Adjustable hangers for most bed types
- High quality pump
- 12 Minute operating cycles
- For patients at a medium to high risk of developing pressure ulcers

2 IN 1 ALTERNATING THERAPY

The Rossmax air mattress range is available in two options, bubble and cell type providing solutions for all four Pressure Ulcer Classifications (Stage 1 – 4).

The bubble type is recommended for stage 1 & 2 where as the cell type is suitable for stage 2-4. The 2 in 1 alternating inflation method inflates half the cells whilst deflating the other half simultaneously.




Anti Decubitus Mattress Bubble Type - AM30




Adjustable hangers to suit any kind of hospital bed frame

Accessories



Order Code: RMAM30



- **2 in 1 alternating therapy**
- **Bubble mattress with durable medical grade PVC**
- **Weight indicator with pressure adjustment settings**
- **Low operating sound**
- Adjustable hangers for most bed types
- Extension flags for stabilisation
- High quality pump
- 6 Minute operating cycles
- For patients at low risk of developing pressure ulcers



rossmax
prevention



Digital Bluetooth Scale 150kg - WF262BT



- Healthstyle App
- Bluetooth® low energy technology
- Body fat assessment
- Body muscle measurement
- BMR (Basal Metabolic Rate) calculation
- BMI (Body Mass Index) calculation
- Visceral fat calculation
- Auto-on/off function
- Graduation 100g
- Maximum weight 150kg

CODE:

RMWF262BT
RMWB101

PRODUCT:

Digital Bluetooth Scale - 150kg
Digital Scale - 180kg

Digital Scale 180kg - WB101



4 Sensor Technology



Auto-On function



Big LCD size



Graduation 100g



Maximum weight 180 kg


- 4 Sensor technology
- Auto-on/off function
- Large LCD display
- Graduation 100g
- Maximum weight 180kg



Unit 4, 68 Parramatta Road, Underwood QLD 4119


1300 778 924 info@constar.com.au www.constar.com.au




WORLD'S FIRST
ARTERY CONDITION CHECK
WITH SpO2 & PULSE RATE
PULSE OXIMETER

ACT

ARTERY CHECK TECHNOLOGY



SB210BT



Z5BT



PARR

PULSE ARRHYTHMIA TECHNOLOGY

early stage
HEART HEALTH
SCREENING
 with ONE single
 measurement



Unit 4, 68 Parramatta Road, Underwood QLD 4119
 1300 778 924 info@constar.com.au www.constar.com.au