



## Air Mattress Bubble

Model: AM30



Adjustable hangers to suit any kind of hospital bed frame

2 to 1 alternating therapy




- Bubble mattress
- 2 to 1 alternating therapy
- Adjustable hangers for most bed types
- Durable medical grade PVC
- Extension flags for stabilization
- High quality pump
- Weight indicator pressure adjustment setting
- 6 minutes cycling time
- Near silent low noise
- For prevention and treatment of patients at low risk of pressure ulcers



Low noise design



Medical grade PVC



2 to 1 alternating therapy



Extension flags for stabilization



Easy to use adjustment



Adjustable hangers to suit any kind of hospital bed frame

### Accessories



High Quality Rubber Tube



fuse



Repairing Kit

Model	Qty per carton	Carton volume
AM 30	4 pcs	0.0408 cbm/ctn

Storage and Transportation Condition  
-18°C to 43°C (0°F to 110°F) 10~90% RH

Distributed by:  
**CONSTAR X PTY LTD**  
Unit 4, 68 Parramatta Road, Underwood, QLD 4119,  
Australia. Ph: 1300 778 924  
[www.constar.com.au](http://www.constar.com.au)