

Air Mattress Cell

Model: AM40




- Cell mattress
- 2 to 1 alternating therapy
- Adjustable hangers for most bed types
- Comfortable Nylon mattress
- High quality pump
- Weight indicator pressure adjustment setting
- 12 minutes cycling time
- Near silent low noise
- For treatment of patients at medium to high risk of pressure ulcers



Low noise design



Comfortable
Comfortable Nylon
mattress



2 to 1 alternating
therapy





Easy to use adjustment



Adjustable hangers to
suit any kind of
hospital bed frame

Adjustable hangers to suit any kind of hospital bed frame

2 to 1 alternating therapy



Accessories



High Quality Rubber Tube



fuse



Spare Cell

Model	Qty per carton	Carton volume
AM 40	2 pcs	0.0511cbm/ctn

Storage and Transportation Condition
-18°C to 43°C (0°F to 110°F) 10~90% RH

Distributed by:

CONSTAR X PTY LTD

Unit 4, 68 Parramatta Road, Underwood, QLD 4119,
Australia. Ph: 1300 778 924

www.constar.com.au