



Infrared Non-Contact Temple Thermometer

Model : HA500 Cutie

Temple Thermometer
Non Contact



Backlight



Storage Case



- 1 second measurement
- Temple Temperature
- Object temperature
- Fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



Patent Nos :
CN: ZL 2011 3 0359484.6
Other patents pending

Clinically Accurate

ASTM E1965-98:2009

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

| Model | Gift box | Qty per carton | Carton volume |
|-------------|----------|----------------|---------------|
| HA500 Cutie | V | 24 pcs | 0.055 cbm/ctn |

Distributed by:
CONSTAR X PTY LTD
 Unit 4, 68 Parramatta Road, Underwood, QLD 4119,
 Australia. Ph: 1300 778 924
 www.constar.com.au

Storage and Transportation Condition
 -25°C ~ 55°C (-13°F ~ 131°F)
 ≤95% RH, 700~1060 hPa