



## Infrared Non-Contact Temple Thermometer

Model: HA500



Temple Thermometer Non Contact



Storage Case



- Measurement distance within 5cm (2 inch)
- 1 second measurement
- Temple Temperature
- Object temperature
- Fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



Model	Qty per carton	Carton volume
HA500	50 pcs	0.072cbm/ctn

Storage and Transportation Condition  
-25°C~55°C(-13°F~131°F)  
≤95% RH, 700~1060 hpa

Patent Nos : ZL 2011 3 0359484.6

### Clinically Accurate

ASTM E1965-98:2009

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

Distributed by:

**CONSTAR X PTY LTD**

Unit 4, 68 Parramatta Road, Underwood, QLD 4119, Australia. Ph: 1300 778 924

[www.constar.com.au](http://www.constar.com.au)