




## Infrared Non-Contact Temple Thermometer

Model: HC700



Patent Nos :  
ZL 2014 3 0054262.7

- Bluetooth low energy Technology
- Telephoto temperature measurement (distant inch)
- 1 second measurement
- Temple temperature
- Object temperature
- Auditory fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



### Clinically Accurate

**ASTM E1965-98:2009**

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.



Model	Qty per carton	Carton volume
HC700	30 pcs	0.060 cbm/ctn

Storage and Transportation Condition  
-25°C~55°C(-13°F~131°F) ≤ 85% RH, 700~1060 hPa

Distributed by:

**CONSTAR X PTY LTD**

Unit 4, 68 Parramatta Road, Underwood, QLD 4119, Australia. Ph: 1300 778 924

www.constar.com.au