

Heavy Duty Wheelchair - Steel Frame

TASK
MEDICAL

The Task Medical heavy duty wheelchair is designed to provide reliable mobility support for patients with higher weight requirements across clinical and care environments. Built with a reinforced double cross brace steel frame and extra-strong adjustable components, it ensures patient comfort while allowing ease of use for caregivers. Its robust construction and 150 kg capacity make it suitable for demanding everyday medical use.

Product Features

- Double cross braces with powder-coated frame
- 8" x 2" PU front tyres
- 24" spoke rear wheels with PU tyres
- Flip-back armrests with desk-length PVC armpad
- Swing-away footrest
- Flame-retardant nylon seat and backrest with sponge padding (extra strong)
- Steel front fork
- Nylon footplate

CODE
TKWHEB



Product Specifications

- | | |
|---------------------------------|-------------------------------|
| ■ Seat Width: 22" | ■ Seat Depth: 45 cm |
| ■ Overall Length: 103 cm | ■ Back Height: 44 cm |
| ■ Overall Width: 78 cm | ■ Floor to Seat Height: 50 cm |
| ■ Overall Height: 91 cm | ■ Weight Capacity: 150 kg |
| ■ Overall Folded Width: 26.5 cm | ■ Front Wheel: 8" |
| ■ Seat Width (cm): 55 cm | ■ Rear Wheel: 24" |

Carton Size

- 88x28x96cm
- G.W 23.0kg
- N.W 20.5kg



www.constar.com.au

