

WHY CHOOSE US?

TASK
MEDICAL

- ✓ Quality You Can Trust
- ✓ Trusted by Healthcare Professionals
- ✓ Comprehensive Product Range
- ✓ Competitive Pricing
- ✓ Reliable Australian Stock Holding
- ✓ Fast Nationwide Delivery
- ✓ Local Customer Support
- ✓ Innovative Healthcare Solutions



Probe covers



Ultrasound gel



Modesty Sheet




Bed Rolls

TASK
MEDICAL



PPE



Biopsy needle



Gynaecology Bed



Instrument Trolley

TASK
MEDICAL

ClearView
Vaginal
Speculum



Safe.
Non-Sterile.
Single-Use S / M / L

Constar X Pty Ltd

1300 778 924

www.constar.com.au

info@constar.com.au



Built to Perform.



Patient Comfort

- Ergonomic design
- Smooth insertion
- Speciality sizing

Clinical Performance

- Clear visual field
- Integrated illumination
- One-handed operation



Safety & Hygiene

- Single-use
- Zero cross contamination
- Latex free

Available in three sizes



Small

Nulliparous patients



Medium

Standard adult use




Large

as clinically indicated




- Clear visualisation with integrated LED illumination designed for safe, efficient examinations.

TASK
MEDICAL



- Smooth, ergonomic design for improved patient comfort.

Instructions for use :



* This is a regulated medical device IFU .

Storage :

- Store in a cool, dry place away from direct sunlight
- Store in original, unopened packaging until point of use
- Do not use after the expiry date printed on the pack



TKSPECDS



TKSPECDM



TKSPEC DL