

Standard Wheelchair - Steel Frame

TASK
MEDICAL

The Task Medical wheelchair is designed to provide reliable mobility support across clinical and care environments. Built with a durable steel frame and practical adjustable components, it ensures patient comfort while allowing ease of use for caregivers. Its robust construction and functional design make it suitable for everyday medical use.

Product Features

- Single cross brace with powder-coated frame
- 8" x 1" PVC front wheels
- 24" spoke rear wheels with pneumatic / PU tyres
- Flip-back armrests with desk-length padded support
- Swing-away footrest with leg support
- Flame-retardant nylon seat and backrest with sponge padding
- Steel front fork for durability
- Plastic footplate

CODE
TKWHE



Product Specifications

- | | |
|--------------------------|---------------------------|
| ■ Seat Width: 18" | ■ Seat Depth: 40 cm |
| ■ Overall Length: 103 cm | ■ Back Height: 41 cm |
| ■ Overall Width: 63 cm | ■ Front to Rear Height |
| ■ Overall Height: 91 cm | 50 cm |
| ■ Overall Folded Width: | ■ Weight Capacity: 120 kg |
| 26.5 cm | ■ Front Wheel: 8" |
| ■ Seat Width (cm): 45 cm | ■ Rear Wheel: 24" |

Carton Size

- 81x27x94cm
- G.W 21.0kg
- N.W 18.5kg



www.constar.com.au

