

Hospital Bed (Electric 3 Function)

TASK
MEDICAL

A dependable three-function electric hospital bed built for clinical environments where patient safety, comfort, and ease of nursing are paramount. Constructed with premium cold-rolled steel and full-length HDPE side rails, it combines practical features with thoughtful safety design for smooth, quiet, reliable everyday use.

Product Features

- Three noise-free motors, remote control operation
- Auto-contour function for patient positioning
- Detachable headboard/footboard, doubles as CPR board
- HDPE side rails with slow-down damping
- Blue liquid angle indicators on side rails
- Infusion holes and drainage hooks
- Height-adjustable IV pole included
- PU sponge mattress included

CODE
TKHB3WG



Product Specifications

- Backrest adjustment: 0–70° ($\pm 5^\circ$)
- Legrest adjustment: 0–40° ($\pm 5^\circ$)
- Height adjustment: 450–750mm (± 10 mm)
- Frame: Cold-rolled steel, electrostatic powder coating
- Castors: Four 125mm, central brake
- Mattress: 1980 x 890 x 80mm



www.constar.com.au

